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Thank you for choosing this upper arm blood pressure monitor from Terraillon. This device allows you to control your blood pressure. It is not intended to be a diagnostic device. Always contact your doctor for any advice, diagnosis or treatment.

CAUTION

- This device is intended for indoor, home use.
- This device is not intended for public use.
- This device is portable, but it is not intended for use during patient transport.
- This device is not suitable for continuous monitoring during medical emergencies or operations.
- This device is intended for no-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than the arm, or for any purpose other than obtaining a blood pressure measurement.
- This device is for adults. Do not use this device on neonates or infants. Do not use it on children unless otherwise instructed by a medical professional.
- Do not use on the women in pregnant, including pre-eclamptic, patients.
- The device is not suitable for use on patients with implanted, electrical devices, such as cardiac pacemakers, defibrillators.
- The effectiveness of this device has not been established for use:
 - on users with common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation,
 - on users with peripheral arterial disease,
 - on users undergoing intravascular therapy, or with arteriovenous (AV) shunt.
 Consult a medical professional before use.
- Do not use this device for diagnosis or treatment of any health problem or disease. Contact your physician if you have or suspect any medical problem. Do not change your medications without the advice of your physician or health care professional.
- If you are taking medication, consult your physician to determine the proper time to measure your blood pressure.

- This device may be used only for the intended use described in this manual, the manufacturer shall have no liability for any incidental, consequential, or special damages caused by misuse or abuse.
- Report any unexpected operation or events to the manufacturer.
- Do not apply the cuff on an arm that has an intravenous drip or a blood transfusion attached.
- Warning: Do not kink, fold, stretch, compress, or otherwise deform the tube during measuring, as the cuff pressure might continuously increase, which could prevent blood flow and result injury.
- Warning: Taking blood pressure measurements too frequently could disrupt blood circulation and cause injuries.
- Warning: Do not apply cuff to areas on patient where skin is delicate or damaged. Check cuff site frequently for irritation.
- Warning: Do not place the cuff on the arm of a person whose arteries or veins are undergoing medical treatment, i.e. intra-vascular access or intra-vascular therapy or an arteriovenous (A-V) shunt, which could disrupt blood circulation and cause injuries.
- Do not place the cuff on the arm on the same side of a mastectomy (especially when lymph nodes have been removed). it is recommended to take measurements on the unaffected side.
- Do not wrap the cuff on the same arm to which another monitoring device is applied. One or both devices could temporarily stop functioning if you try to use them at the same time.
- Please check that the operation of the device do not result in prolonged impairment of patient blood circulation.
- Warning: On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, loosen and remove the cuff immediately. Prolonged high pressure applied to the arm (cuff pressure >300 mmHg or constant pressure >15 mmHg for more than 3 minutes) might lead to bruising and discolored skin.
- Warning: Do not use this device with high-frequency (HF) surgical equipment at the same time.
- Warning: This device is not AP/APG equipment. Do not use the device where flammable anesthetic are present, or in environments mixture with air of with oxygen or nitrous oxide.
- The device contains sensitive electronic components. To avoid measurement errors, avoid taking blood pressure measurements near a strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.

- Wireless communication equipment, such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies may cause interference that may affect the accuracy of measurements. A minimum distance of 1 foot (30 cm) should be kept from such devices during a measurement.
- You can use this device to take your own measurement, no third-party operator is required.
- Please use the device under the environment which is provided in the user manual. Otherwise, the performance and lifetime of the device will be impacted and reduced.
- The device may require up to 30 minutes to warm up / cool down from the minimum / maximum storage temperature before it is ready for use.
- Warning: Excessive cuff tube lengths could cause strangulation if you don't manage them properly.
- Warning: Do not touch output of the batteries/adaptor and the user simultaneously.
- Adaptor is specified as a part of ME EQUIPMENT.
- Warning: The power cord is considered the disconnect device for isolating this equipment from supply mains. Do not position the equipment so that it is difficult to reach or disconnect.
- The blood pressure monitor, its adapter, and the cuff are suitable for use within the patient environment.
- Warning: Do not use this device if you are allergic to polyester, nylon, or plastic.
- Warning: Only use accessories approved by manufacturer. Using unapproved accessories might cause damage to the unit and injure users.
- Warning: If you experience discomfort during a measurement, such as pain in the arm or other complaints, press the Power button immediately to release the air from the cuff.
- No calibration is required within two years of reliable service.
- Do not attempt to repair the unit yourself if it malfunctions. Only have repairs carried out by authorized service centers.
- At the request of authorized service personnel, circuit diagrams, component part lists, descriptions, and calibration procedures will be made available by the manufacturer or distributor.
- It is recommended that the performance should be checked after repair, maintenance, and every two years of use, by retesting the requirements in limits of the error of the cuff pressure indication and air leakage (testing at least at 50 mmHg and 200 mmHg).

- This device may be used only for the intended use described in this manual, the manufacturer shall have no liability for any incidental, consequential, or special damages caused by misuse or abuse.
- Report any unexpected operation or events to the manufacturer.
- Do not apply the cuff on an arm that has an intravenous drip or a blood transfusion attached.
- Warning: Do not kink, fold, stretch, compress, or otherwise deform the tube during measuring, as the cuff pressure might continuously increase, which could prevent blood flow and result injury.
- Warning: Taking blood pressure measurements too frequently could disrupt blood circulation and cause injuries.
- Warning: Do not apply cuff to areas on patient where skin is delicate or damaged. Check cuff site frequently for irritation.
- Warning: Do not place the cuff on the arm of a person whose arteries or veins are undergoing medical treatment, i.e. intra-vascular access or intra-vascular therapy or an arteriovenous (A-V) shunt, which could disrupt blood circulation and cause injuries.
- Do not place the cuff on the arm on the same side of a mastectomy (especially when lymph nodes have been removed). It is recommended to take measurements on the unaffected side.
- Do not wrap the cuff on the same arm to which another monitoring device is applied. One or both devices could temporarily stop functioning if you try to use them at the same time.
- Please check that the operation of the device do not result in prolonged impairment of patient blood circulation.
- Warning: On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, loosen and remove the cuff immediately. Prolonged high pressure applied to the arm (cuff pressure >300 mmHg or constant pressure >15 mmHg for more than 3 minutes) might lead to bruising and discolored skin.
- Warning: Do not use this device with high-frequency (HF) surgical equipment at the same time.
- Warning: This device is not AP/AG equipment. Do not use the device where flammable anesthetic are present, or in environments mixture with air of with oxygen or nitrous oxide.
- The device contains sensitive electronic components. To avoid measurement errors, avoid taking blood pressure measurements near a strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.

Please check that the operation of the device does not cause prolonged impairment of the patient's blood flow.

In the rare event of a fault that causes the cuff to remain fully inflated during measurement, loosen, and remove the cuff immediately. Prolonged high pressure applied to the arm (cuff pressure >300 mmHg or constant pressure >15 mmHg for more than 3 minutes) may cause bruising and skin discoloration.

Performance may be affected by extremes of temperature, humidity and altitude.

The unit may need 30 minutes to warm up or cool down from the minimum or maximum storage temperature before being ready for use.

Make sure the rubber cuff tube is not crushed, stretched, or bent during storage.

CONTRAINDICATIONS


1. The device should not be used by any person who may be suspected of, or is pregnant .
2. The device is not suitable for use on patients with implanted, electrical devices, such as cardiac pacemakers, defibrillators.

CAUTION:


- The device is not suitable for use on neonatal patients, pregnant women, patients with implanted, electrical devices, patients with pre-eclampsia, premature ventricular beats, atrial fibrillation, peripheral, arterial disease and patients undergoing intravascular therapy or arterio-venous shunt or people who received a mastectomy. Please consult your doctor prior to using the unit if you suffer from illnesses.
- The device is not suitable for measuring the blood pressure of children. Ask your doctor before using it on older children.
- The device is not intended for patient transport outside a healthcare facility.
- The device is not intended for public use.
- When using this device, please pay attention to the following situation which may interrupt blood flow and influence blood circulation of the patient, thus cause harmful injury to the patient: connection tubing kinking too frequent and consecutive multiple measurements; the application of the cuff and its pressurization on any wrist where intravascular access or therapy, or an arterio-venous (A-V) shunt, is present; inflating the cuff on the side of a mastectomy.

- Warning: Do not apply the cuff over a wound; otherwise it can cause further injury.
- The maximum temperature that the applied part can be achieved is 42.5°C while the environmental temperature is 40°C.
- Before use, make sure the device functions safely and is in proper working condition. Check the device; do not use the device if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.
- Do not wash the cuff in a washing machine or dishwasher!
- At least 30 min required for ME equipment to warm from the minimum storage temperature between uses until it is ready for intended use. At least 30 min required for ME equipment to cool from the maximum storage temperature between uses until it is ready for intended use.
- This equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS;
- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance d away from the equipment. The distance d is calculated by the MANUFACTURER from the 80MHz to 5.8 GHz column of Table 4 and Table 9 of IEC 60601-1-2:2014, as appropriate.
- Please use the device under the environment which was provided in the user manual. Otherwise, the performance and lifetime of the device will be impacted and reduced.

MONITOR COMPONENTS

LCD DISPLAY SIGNAL


SYMBOL	DESCRIPTION	EXPLANATION
SYS	Systolic blood pressure	High pressure result
DIA	Diastolic blood pressure	Low pressure result
Pul/min	Pulse per minute	Heartbeats per minute, BPM
▼	Deflating	CUFF air is exhausting of deflating
MEM	Memory	The displayed measurement values is from the memory
kPa	kPa	Measurement Unit of the blood pressure
mmHg	mmHg	Measurement Unit of the blood pressure
Lo +	Low battery	Batteries are low and need to be replaced
Irregular heartbeat	Irregular heartbeat	Irregular heartbeat detection
Grade	Grade	The grade of the blood pressure
Time	Current Time	Year/Month/Day, Hour/Minute
A	User A	Start measurement for User A
B	User B	Start measurement for User B
Heartbeat	Heartbeat	Heartbeat detection during the measurement
AVG	The average value	The average value of the latest three records

BEFORE YOU START

Installing and Replacing the Batteries

- 4×AAA batteries (Not Included)
- Open the battery cover
- Install the batteries by matching the correct polarity, as shown
- Replace the cover

Replace the batteries whenever the below happen

- The  shows
- The display dims
- The display does not light up

CAUTION

Remove batteries if the device is not likely to be used for some time. The old batteries are harmful to the environment; do not dispose with other daily trash. Remove the old batteries from the device and follow your local recycling guidelines. Do not dispose of batteries in fire. Batteries may explode or leak.

Setting Date, Time and Measurement Unit


It is important to set the clock before using your blood pressure monitor, so that a time stamp can be assigned to each record that is stored in the memory.

1. When the monitor is off, hold pressing "SET" button for 3 seconds to enter the mode for year setting. Or when the monitor is off, press "SET" button shortly, it will display the time. Then hold pressing "SET" button to enter the mode for year setting. Note: "MEM" button allows you to change the data, "SET" button allows you to confirm the data.
2. Press "MEM" button to change the year [Yr].
3. When you get the right year, press "SET" button to set down and turn to next step.
4. Repeat steps 2 and 3 to set the month [M] and day [D].
5. Repeat steps 2 and 3 to set the hour and minute.
6. Repeat steps 2 and 3 to set the unit.
7. After the unit is set, the LCD will display "dOnE" first, then display all the settings you have done and then it will turn off.

MEASUREMENT

APPLYING THE CUFF

1. Remove all jewelry, such as watches and bracelets, from your left arm.
Note: If your doctor has diagnosed you with poor circulation in your left arm, use your right arm.
2. Roll or roll up your sleeve to expose the skin. Make sure your sleeve is not too tight.
3. Hold your arm with your palm up and apply the cuff to your upper arm, then position the tube off-center towards your inner arm, in line with your little finger or position the artery marker on the main artery (inside your arm). Note: Locate the main artery by pressing with 2 fingers about 2 cm above the elbow crease on the inside of the left arm. Identify where the pulse is strongest. This is your main artery.
4. The cuff should fit snugly but not too tightly. You should be able to insert a finger between the cuff and your arm.
5. Sit comfortably with your arm on a flat surface. Place your elbow on a table so that the cuff is level with your heart. Turn your palm upwards. Sit upright in a chair and take 5-6 deep breaths.
6. Useful tips for users, especially for users with high blood pressure:
 - Rest for 5 minutes before the first measurement
 - Wait at least 3 minutes between measurements. (This allows your blood circulation to recover). Take the measurement in a quiet room.
 - The user should relax as much as possible and not move or talk during the measurement procedure.
 - The cuff should be positioned at the same level as the right atrium of the heart.
 - Sit comfortably. Do not cross your legs and keep your feet flat on the floor.


- Keep your back against the back of the chair.
- For a meaningful comparison, try measuring under similar conditions. For example, take daily measurements at about the same time, on the same arm, or as instructed by the doctor.

Start the Measurement

- Before you start the measurement, switch the User button to select the user between User A and User B. Switch to right to select User A, switch to left to select User B. When the monitor is off, press the "START/STOP" button to turn on the monitor, and it will finish the whole measurement. And save the measurement data for the desired user. (Take User A for example.)

- (1) LCD display
- (2) Adjust to zero
- (3) Inflating and measuring
- (4) Display and save the measurement result

- Press "START/STOP" button to power off, otherwise it will turn off within 1 minute.

Tips: Maximum 250 records are both for User A and User B.

DATA MANAGEMENT

Recall the Records

- When the monitor is off, please press "MEM" button, it will display the latest record first when the records are less than three groups. When there are more than three groups, it will display the average value of the latest three records first.
- Press "MEM" button or "SET" button to get the record you want. Press and hold "MEM" button to look over ten groups of the historical records quickly.

The date and time of the record will be shown alternately.

The current No. is No 2

The corresponding date
is January 2nd

The corresponding time
is 10:08

- If you want to look over another user's data, switch the User button to select the desired user. Then you can look over its historical records.

Delete the Records

If you did not get the correct measurement, you can delete all results for the selected user by following steps below.

- Hold pressing "SET" button for 3 seconds when the monitor is in the memory recall mode, the flash display "dEL ALL" will show.
- Press "SET" button to confirm deleting and the monitor will turn off.
- If you don't want to delete the records, press "START/STOP" button to escape.

INFORMATION FOR USER

Tips for Measurement

Measurements may be inaccurate if taken in the following circumstances:

- Within 1 hour after dinner or drinking
- Immediate measurement after tea, coffee, smoking
- Within 20 minutes after taking a bath
- When talking or moving your fingers
- In a very cold environment
- When you want to discharge urine

Maintenance


In order to get the best performance, please follow the instructions below:

- Put in a dry place and avoid the sunshine
- Avoid touching water, clean it with a dry cloth in case
- Avoid intense shaking and collisions
- Avoid dusty and unstable temperature environment
- Using wet cloths to remove dirt
- Do not attempt to clean the reusable cuff with water and never immerse the cuff in water

ABOUT BLOOD PRESSURE

What are systolic pressure and diastolic pressure?

When ventricles contract and pump blood out of the heart, the blood pressure reaches its maximum value in the cycle, which is called systolic pressure. When the ventricles relax, the blood pressure reaches its minimum value in the cycle, which is called diastolic pressure.

What is the standard blood pressure classification?

The blood pressure classification published by World Health Organization (WHO) and International Society of Hypertension (ISH) in 1999 is as follows:

CAUTION

Only a physician can tell your normal BP range. Please contact a physician if your measuring result falls out of the range. Please note that only a physician can tell whether your blood pressure value has reached a dangerous point.

≥180	severe						
160 ~ 179	medium						
140 ~ 159	moderate						
130 ~ 139	high normal						
120 ~ 129	normal						
< 120	optimum						
SYSTOLIC mmHG		optimum < 80	normal 80~84	high normal 85~89	moderate 90~99	medium 100~109	severe ≥110
DIASTOLIC mmHG							

Irregular Heartbeat Detector

An irregular heartbeat is detected when a heartbeat rhythm varies while the device is measuring systolic pressure and diastolic pressure. During each measurement, blood pressure monitor will keep a record of all the pulse intervals and calculate the average value of them. If there are two or more pulse intervals, the difference between each interval and the average is more than the average value of $\pm 25\%$, or there are four or more pulse intervals, the difference between each interval and the average is more than the average value of $\pm 15\%$, then the irregular heartbeat symbol will appear on the display with the measurement result.

CAUTION

The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

Why does my blood pressure fluctuate throughout the day?

1. Individual blood pressure varies multiple times everyday. It is also affected by the way you tie your cuff and your measurement position, so please take the measurement under the same conditions.
2. If the person takes medicine, the pressure will vary more.
3. Wait at least 3 minutes for another measurement.

Why do I get a different blood pressure at home compared to the hospital?

The blood pressure is different even throughout the day due to weather, emotion, exercise etc. Also, there is the “white coat” effect, which means blood pressure usually increases in clinical settings.

What you need to pay attention to when you measure your blood pressure at home:

- If the cuff is tied properly on the upper arm.
- If the cuff is too tight or too loose.
- If you feel anxious.
- Taking 2-3 deep breaths before beginning will be better for measuring.
- Advice: Relax yourself for 4-5 minutes until you calm down.

Is the result the same if measuring on the right arm?

It is ok for both arms, but there will be some different results for different people. We suggest you measure the same arm every time.

TROUBLESHOOTING

PROBLEM	SYMPTOM	CHECK THIS	REMEDY
No power	Display will not light up	Batteries are exhausted	Replace with new batteries
		Batteries are inserted incorrectly	Insert the batteries correctly
Low batteries	Display is dim or show 	Batteries are low	Replace with new batteries
Error message	E 1 shows	The cuff is not secure	Refasten the cuff and then measure again
	E 3 shows	The pressure of the cuff is excess	Relax for a moment and then measure again
	E10 ou E11 shows	The monitor detected motion, talking or the pulse is too poor while measuring	Relax for a moment and then measure again
	E20 shows	The measurement process does not detect the pulse signal	Loosen the clothing on the arm and then measure again
	E21 shows	The treatment of the measurement failed	Relax for a moment and then measure again
	EExx,shows on the display	A calibration error occurred (XX can be some numeric characters,such as 01, 02 and so on)	Retake the measurement. If the problem persists, contact the retailer or our customer service department for further assistance.Refer to the warranty for contact information and return instructions
Warning message	“out ” shows	Out of measurement range	4.Relax for a moment. Refasten the cuff and then measure again. If the problem persists, contact your physician.

WARRANTY & ENVIRONMENT PROTECTION

This product is guaranteed for 2 years against material and manufacturing defects. During this period, product will be replaced free of charge (proof of purchase must be presented if the product is under warranty). This warranty does not cover damage resulting from accidents, improper use or negligence. If you have a complaint please contact the store where you purchased the product first.

Waste electrical products should not be disposed of with household waste.

Please recycle where facilities exist.

Check with your local authority or retailer for recycling advice.

SAFETY INFORMATION

The signs below might be in the user manual, labeling or other component. They are the requirement of standard and using.

	Refer to instruction manual/booklet . To signify that the instruction manual/ booklet must be read.		Type BF applied part
CE 0123	CE marking indicates that a product has been assessed by the manufacturer and deemed to meet EU safety, health and environmental protection requirements. It is required for products manufactured anywhere in the world that are then marketed in the EU.		The symbol indicates that the product should not be discarded as unsorted waste but must be sent to separate collection facilities for recovery and recycling.
	The name and the address of the authorised representative in Switzerland		Authorized representative in the European Community/European Union
	Symbol for "SERIAL NUMBER"		Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences
	Symbol for "DIRECT CURRENT"		
	Date and Country of manufacture		
	The packaging is recyclable		
	Symbol for "MANUFACTURER"		

SPECIFICATIONS

Power supply	Battery powered mode: 6VDC 4×AAA batteries (Not Included). AC adaptor powered mode: 6V  1A (Please only use the recommended AC adaptor model). (Not Included)
Display mode	Digital LCD V.A.60mm×92mm
Measurement mode	Oscillographic testing mode
Measurement range	Rated cuff pressure: 0kPa - 40kPa (0mmHg~300mmHg). Measurement pressure: 5.3kPa-30.7kPa (40mmHg-230mmHg) pulse value: (40-199) beat/minute
Accuracy	Pressure: 5°C-40°C within±0.4kPa(3mmHg) pulse value:±5%
Normal working condition	A temperature range of: +5°C to +40°C A relative humidity range of 15% to 90%, non-condensing, but not requiring a water vapour partial pressure greater than 50 hPa An atmospheric pressure range of : 700 hPa to 1060 hPa

Storage & transportation condition	Temperature:-20°C to +60°C A relative humidity range of ≤ 93%, non-condensing, at a water vapour pressure up to 50hPa
Measurement perimeter of the upper arm	About 22cm~42cm
Net Weight	Approx.250g (Excluding the dry cells)
External dimensions	Approx.140mm×130mm×49.7mm
Attachment	user manual
Mode of operation	Continuous operation
Degree of protection	Type BF applied part
Protection against ingress of water	IP 21:It means the device could be protected against solid foreign objects of 12,5mm ø and greater, and against vertically falling water drops.
Software Version	A01
Device Classification	Battery Powered Mode: Internally Powered ME Equipment AC Adaptor Powered Mode: Class II ME Equipment

WARNING: No modification of this equipment is allowed.

MEASUREMENT PRINCIPLE

This product uses the Oscillometric Measuring method to detect blood pressure.Before every measurement, the unit establishes a "zero pressure" equivalent to the air pressure. Then it starts inflating the arm cuff, meanwhile, the unit detects pressure oscillations generated by beat-to-beat pulsatile, which is used to determine the systolic and diastolic pressure, and also pulse rate.

INDICATIONS FOR USE

The Blood Pressure Monitor is digital monitors intended for use in measuring blood pressure and heartbeat rate with arm circumference ranging from 22cm to 32 cm (about 8¾"-12½") or 22cm to 42cm(about 8¾"-16½"). It is intended for adult indoor use only.

CONTRAINDICATIONS

- 1.The device should not be used by any person who may be suspected of,or is pregnant .
- 2.The device is not suitable for use on patients with implanted,electrical devices, such as cardiac pacemakers, defibrillators.

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Terrailon.

EMC GUIDANCE

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments.

Warning: Don't be near the active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment TMB-1490-S including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Technical description:

1. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected service life.
2. Guidance and manufacturer's declaration-electromagnetic emissions and Immunity.

Table 1

Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations / flicker emissions IEC 61000-3-3	Comply

Table 2

Guidance and manufacturer's declaration – electromagnetic Immunity		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV signal input/output 100 kHz repetition frequency	±2 kV for power supply lines Not Applicable 100 kHz repetition frequency
Surge IEC61000-4-5	±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV, ±2 kV common mode	±0.5 kV, ±1 kV differential mode Not Applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT; 1 cycle and 70% UT; 25/30 cycles; Single phase: at 0°. 0% UT; 250 / 300 cycle	0% UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT; 1 cycle and 70% UT; 25/30 cycles; Single phase: at 0°. 0% UT; 250 / 300 cycle
Power frequency magnetic field IEC 61000-4-8	30 A/m 50 Hz / 60 Hz	30 A/m 50 Hz / 60 Hz
Conducted RF IEC61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz
NOTE U _i is the a.c. mains voltage prior to application of the test level.		

Table 3

Guidance and manufacturer's declaration - electromagnetic Immunity							
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
	385	380-390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27
	450	430-470	GMRS 460, FRS 460	FM $\pm$ 5kHz deviation 1kHz sine	2	0.3	28
	710	704-787	LTE Band 13, 17	Pulse modulation 217Hz	0.2	0.3	9
	745						
	780						
	810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	2	0.3	28
	870						
	930						
	1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217Hz	2	0.3	28
	1845						
	1970						
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
	5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
	5500						
	5785						

Terraillon SAS France & Headquarters

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Terraillon Asia Pacific Ltd

4/F, Eastern Centre
1065 King's Road
Quarry Bay - Hong Kong
Tel: + 852 (0)2960 7200
customerservice@terraillon-asia.com

Terraillon Corp USA

contact@terraillon.fr

Terraillon®

www.terraillon.com