TALIKA

PARIS

TALIKA
PARIS

HAIR FORCE LED CAP

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EN

A pioneer in electro-beauty for over 20 years, Talika presents HAIR FORCE LED CAP, the ultimate ANTI-HAIR LOSS and FORTIFYING solution for healthy, stronger, denser and longer hair. Lightweight, ergonomic and stylish, this cap you can take everywhere with you revitalizes the scalp and stimulates growth thanks to its 120 red LEDs at 625nm.

STIMULATING HAIR GROWTH HAS NEVER BEEN SO EASY !

- Reduces hair loss
- Reduces scalp micro-inflammations
- Boosts microcirculation
- Stimulates hair growth
- Increases hair density

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HAIR FORCE LED CAP

Parameters

LED WAVELENGTHS	625+/-5 NM (RED)
Max. output optical power density	10 to 20 mW/cm ²
LEDs	120 Red LEDs
Operating temperature	10°C to 40°C 50°F to 104°F
Storing and transport temperature	-5°C to 55°C 23°F to 131°F
Operating humidity	30% to 95% without condensation
Storage and Transport Humidity	20% to 95% without condensation
Charging supply	5V 1A
Charging port	USB-C
Charging cable	USB-C to USB-A
Charging adapter power output	1A
Battery	3.7V, 1800 mAh
Battery type	Lithium-ion
Charging time	240 minutes / 4 hours
Recommended use	Once in 24 hours
Treatment time	10 minutes
Dimensions of the device	295 x 170 x 122 mm
Weight of the device	245 g / 8.64 oz

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RED LIGHT - 625 NM:

This wavelength is known for its soothing, regenerating and stimulating effects.
By increasing collagen and elastin synthesis, it improves the quality of the scalp to promote and limit the effects of aging.
It also reduces local micro-inflammations.

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EN Preparation before using

Remove all products from your hair and scalp and making sure your hair and scalp are completely dry. Make sure the device is unplugged and place it on your head. Adjust the cap to your size if necessary.

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EN Use the device once a day for a 3-month cure. It is not recommended to use the device more than once a day.

Fully charge your device for 2 hours before using it for the first time. Only use a 5V adapter or the USB port on your computer (avoid fast chargers as they can damage the battery).

The LED on the device will start to flash slowly to indicate charging. When the device is fully charged, the LED indicator will go out.

Before charging the device, make sure that the components are dry. Failure to do so may result in an electric shock, short-circuit or fire.

Note: always unplug the device before using it, as it will not work while it is charging.

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EN We highly recommend users perform a light sensitivity test before their first use of the device.

1. Place the device on the inner forearm.
2. Turn on the device to start the therapy.
3. Wait till the treatment is done. The treatment will take 10 minutes to complete.
4. If your skin feels hot or uncomfortable during the treatment directly stop! You might be sensitive to light.
5. The device will automatically turn off after the treatment is done. Remove the device and wait for at least 24 hours.

Check your skin regularly over those 24 hours, if your skin appears red and/or itchy the device is unsuitable for you. In the case of light sensitivity please contact customer support. (see last page of manual).

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EN If you experience any adverse effects, such as discomfort, skin discolouration, skin irritation, headaches or other problems related to the use of this device, stop using it immediately and consult your Healthcare provider.

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EN - Immediately after treatment skin may appear red. Skin should return to normal within 1-2 hours.
- Temporary headaches immediately following the light treatment.

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EN

- The lifetime of the Lithium battery can be maximized by allowing it to fully charge and then fully discharge during regular use.
- Store the device at room temperature in dry conditions.
- Avoid storing the device in dusty conditions, in direct sunlight or in damp and humid conditions.
- Before each use, with the device switched off, clean the Device following the cleaning instructions.
- This device contains no serviceable parts.
- If your device malfunctions in any way, please contact Talika customer service. *(see last page of manual)*.

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EN

- DO NOT dispose of the device with household waste. The Device should be recycled.
- The Device contains a Lithium-ion battery. For protection of the environment, please safely dispose of the device following your local guidelines.

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CONTRAINDICATIONS & WARNINGS

EN

- DO NOT use the device if you have any active cancer or cancerous lesion.
- DO NOT use the device if you are pregnant, planning on becoming pregnant, or breastfeeding.
- DO NOT use the device on people under the age of 18.
- DO NOT use the device over areas with phlebitis, thrombophlebitis, and varicose veins.
- DO NOT use the device if you have a known allergic reaction to **Plastics or Neoprene**. The device contains Plastics, Fabric and Metals that come into contact with the skin during use. When in doubt of allergies first contact your Healthcare Provider before using the device.
- DO NOT use the device if you suffer from any photosensitive disorder (sensitization to light).
- DO NOT use the device if you are taking any medication or drugs that can cause photosensitivity. If you are not sure about your medication, contact your Healthcare Provider.
- DO NOT use the device if you suffer from light-induced headaches.
- DO NOT use the device if you have a medical history of seizures or epilepsy that was triggered by light.

Note: We highly suggest taking caution when taking any (skin related) damaging medications. If in doubt, consult your Healthcare Provider.

Note: If you are under the care of a Healthcare Provider, consult your Healthcare Provider before using this device.

- STOP using the device and consult your Healthcare Provider if you experience an adverse reaction from using the device.
- STOP using the device if you experience symptoms of hypotension (dizziness or light-headedness) while using the device.
- STOP using the device if it is overheating. Turn off the device immediately.
- STOP your treatment if it feels hot or your skin feels uncomfortably warm.
- DO NOT machine wash the device.
- DO NOT let children play with the device. Cleaning and charging the device shall not be done by children.
- DO NOT use the device to treat rosacea, moles, warts, open sores, cancerous lesions, or any skin condition.
- DO NOT shine the light from the LEDs directly into your eyes or someone else's eyes as this may cause temporary loss of sight.
- DO NOT place the device over any open wounds.
- DO NOT sit, lay down or step on the device. This will break the device.
- DO NOT drop or insert any object into any openings of the device.
- DO NOT use the device if it is not working properly, or if it has been dropped or damaged. For warranty and repairs contact customer support.
- DO NOT use or store the device near hot surfaces.
- DO NOT charge the device when it is close to or might get in contact with water.
- DO NOT submerge the device in water.
- DO NOT store the device in extreme heat or cold.
- DO NOT place the device in direct sunlight.
- DO NOT attempt to modify the device in any way. Modifying or opening the device will void your warranty. The device contains no serviceable parts. This device contains batteries that are non-replaceable.
- DO NOT perform maintenance or take apart the product. In case of malfunction or damage, contact customer support.
- DO NOT use the device if damaged or when the charging cable is damaged.
- DO NOT use the device with non-manufacturer-approved accessories.
- DO NOT use solvents or abrasive cleaners on the device.
- DO NOT use the device in any other way than specified in this manual.
- DO NOT use the device if you have had a recent surgery or other surgical procedures.
- DO NOT use the device over a bacterial infection as this may result in the worsening of symptoms.
- DO NOT use the device on painful areas. Consult your Healthcare Provider when you experience painful area's during use.
- DO NOT use the device on, near, or by children or those with reduced physical and mental abilities. Close supervision is necessary when this device is used, cleaned, or stored near children or those with reduced physical or mental abilities.
- DO NOT use the device with any other USB-C cable than provided by Talika.

Note: The device must only be supplied with a low voltage of 5V DC. If any part of the main adapter or cable is damaged, stop using the device and contact Talika customer service (<mailto:contact@talika.com>).

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EN This device is not intended to diagnose, cure or prevent any disease. Talika makes no claims, representations, or warranties regarding the ability of this device to cure any physical, skin, or mental condition. A qualified Healthcare Provider should be consulted regarding any condition requiring medical attention.

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EN The device's lightweight design makes it a perfect travel companion. Here are some travel tips for you:
- If you are travelling to a different region of where you bought the Device, please assure that the device is only supplied with a 5V DC Low Voltage.
- Confirm with the steward that the Device can be carried and/or used on the airplane prior to flight.

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EN. Scan the QR code for more info & activate the warranty.

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