essentiel 

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Thank you for purchasing an Essentiel b brand product. We pay special attention to the **RELIABILITY, EASE OF USE** and **DESIGN** of our products. We hope that you will be fully satisfied with this body fat scale.

## your product


### Contents of the box

- 1 body fat scale
- 2 x 1.5 V AAA batteries
- 1 user manual

### Technical features


- Maximum weight: 150 kg
- Graduation: 100 g / 1 kcal
- LCD screen: 74 x 40 mm
- Height of digits: 30 mm
- Unit of measurement: kg/lb
- 22 mm glass platform
- Uses 2 x 1.5 V AAA batteries
- Automatic switch-on by pressing with your foot or by pressing the settings button
- Auto-off
- Muscle mass
- Fat mass
- Water mass
- Bone mass
- Overload indicator
- Used battery indicator
- 4 anti-slip feet
- 12 user memories
- Dimensions (L x H x W): 302 x 22 x 302 mm

### Appliance description

1	LCD screen
2	Decrease button
3	Settings button
4	Increase button
5	Sensors x4

### Display screen description

### Button functions

**SET** : Settings button. Turns on the device and scrolls through the settings.

▲ : Increases the value by one unit with each press. Pressing and holding increases the values quickly.


▼ : Decreases the value by one unit with each press. Pressing and holding decreases the values quickly.

### User data settings


User data	Setting values
Memory number	12 memories
Sex	Male  , Female 
Height	from 1 metre to 2.50 m (100 to 250 cm)
Age	from 10 to 100 years

## before initial use

### Choice of measurement unit

You can change your weight's unit of measurement. This means you can choose between **kilograms (kg)** or **pounds (lb or st-lb)**<sup>1</sup>. Turn on the device, then press the button at the rear of the device (see diagram opposite); the screen will display the chosen unit. Press again to select a different unit.


Please note:

- When you select kilograms, the default height unit will be centimetres (cm).
- When you select pounds (lb or st-lb), the default height unit will be feet/inches<sup>2</sup>.

### Battery installation/replacement

1. Remove the battery compartment cover located under the body fat scale.
2. Remove used batteries.
3. Replace them with two new 1.5 V AAA batteries, observing the correct polarity.

- Use only the type of battery specified in this manual.
- Replace the batteries when they are worn out or leaking.
- Remove the batteries if the device is not used for a long time.

1 1 English pound = approx. 453 grams. 1 stone = 14 pounds = 6.350 kilograms  
2 1 foot = approx. 30 cm or 12 inches. 1 inch = 2.54 cm.

## useful information

### 1) Fat mass

Fat mass is the percentage of total weight that evaluates the proportion of fat in the body. The calculation of fat mass applies only to adults between 18 and 65 years of age and not to very young or very old people, pregnant or breastfeeding women, very muscular people, high-level athletes, very tall or very short people, and people with limb amputations. The ideal fat mass varies according to age and sex. For all people outside the normal range, a check-up can be carried out with a doctor, nutritionist or dietician.

#### Why is it useful to control the percentage of fat in the body?

The total weight determines whether a person is obese or not, but the change in weight alone does not indicate whether it is a change in fat mass or muscle mass. To control your weight, it is best to maintain your muscle mass while losing fat mass. Controlling the fat percentage is therefore an important step towards optimal weight management and good health.

You can refer to the following table. Every individual is different; this table is provided for information purposes only.

#### Fat mass table

Age	Woman				
≤ 30 years	≤16%	16.1 - 20.5%	20.6 - 25%	25.1 - 30.5%	30.5%
> 30 years	≤ 20%	20.1 - 25%	25.1 - 30%	30.1 - 35.5%	35.5%
On-screen display	Very sufficient	Insufficient	Normal	Above normal	Much higher than normal
	i--	i-	i	i+	i++

Age	Man				
≤ 30 years	≤11%	11.1 - 15.5%	15.6 - 20%	20.1 - 24.5%	24.5%
> 30 years	≤15%	15.1 - 19.5%	19.6 - 24%	24.1 - 28.5%	28.5%
On-screen display	Very sufficient	Insufficient	Normal	Above normal	Much higher than normal
	i--	i-	i	i+	i++

Note: The fat percentage calculated by the body fat scale is only an approximation of your actual fat mass. There are precise medical methods that calculate fat mass. Seek medical advice.

## 2) Water mass

Water mass is the percentage of total weight that evaluates the proportion of liquid in the body. The body of a healthy individual consists of about 50-70% water. Water plays an essential role in many functions of the human body and is present in cells, tissues and organs. The ideal water mass varies according to age and sex.

### Why is it useful to control the percentage of water in the body?

A good water mass level indicates the proper functioning of the body and reduces the risk of developing health problems. Water mass varies according to your hydration level (the amount of water you have drunk or have lost in perspiration). For greater accuracy, it is better to measure your water mass in the early evening, before dinner, when your hydration level is most stable. It is not advisable to measure your water mass after getting up in the morning, exercising or bathing because of dehydration and the poor distribution of liquid in your body. In this case, wait until the liquid is more evenly distributed under the effect of gravity.

Factors that influence water mass are sleep, physical activity, eating and drinking.

### Water mass table

The water mass rate (in %) is usually in the following ranges:

Woman			
Age	Poor	Good	Very good
< 40 years	45%	45~60%	60%
40-60 years	43.5%	43.5~58.5%	58.5%
> 60 years	42%	42~57%	57%

Man			
Age	Poor	Good	Very good
< 40 years	50%	50~65%	65%
40-60 years	48.5%	48.5~63.5%	63.5%
> 60 years	47%	47~62%	62%

## 3) Muscle mass

Muscle mass is the percentage of total weight that evaluates the proportion of muscle in the body. Muscles transform food into mechanical energy and maintain body temperature. They are the foundation of our physical condition.

### Why is it useful to control the percentage of muscle in the body?

Controlling muscle mass is particularly important for people with weight problems. The body needs muscles to eliminate excess energy reserves, and effective and long-lasting weight loss must be combined with muscle building.

## 4) Bone mass

Bone mass is the proportion of weight made up of skeletal bones. The building of bone capital begins in childhood and adolescence, periods during which it is important to have a satisfactory calcium intake. Then, having reached an optimum level around the age of 25, this bone capital remains constant for about ten years, before decreasing.

In men, bone loss is then linear and of low intensity, of the order of 0.5% per year. But in women, it accelerates at menopause due to oestrogen deficiency to reach 3 to 5% for 2 to 3 years and then 1 to 2% during the next 10 years.


Osteoporosis is accelerated bone loss. Treatments exist, but the best way to fight this is prevention with a varied and balanced diet.

### User tips

- In the morning, wait about 15 minutes after getting up to take the reading to allow the water to be properly distributed in the body.
- Wait a few hours to take the reading after a physical effort.
- Stand upright and do not move during the reading.
- You should always be barefoot when taking the reading; however, the soles of the feet may be slightly damp. Please note that measurements taken with completely dry feet do give satisfactory results.
- Please note that the results may not be relevant in the following cases: for high-level athletes, people with fever, on dialysis, with symptoms of oedema or osteoporosis, people undergoing cardiovascular medical treatment or taking vasodilator or vasoconstrictor drugs.

## using the scale

### Weighing only

1. Place the scale on a flat and stable surface for safety and accuracy.
2. Gently step onto the platform and stand still until the numbers stabilise.
3. The numbers flash and then remain fixed. Step off the scale. It then turns off automatically.

## using the body fat scale


### Definition of settings

Before the first use, you will need to set up your personal data. To do so, proceed as follows:

#### 1. Select a memory number

You can configure personal data for up to 12 users.


- Turn on the device by pressing the **SET** button to enter the configuration menu. The memory number will blink.
- Press the ▲ or ▼ button to select a memory number.

- Confirm your choice by pressing the **SET** button, you will move to the next setting.

#### 2. Select sex


- Press the ▲ or ▼ button to select: Male  or Female .

- Confirm your choice by pressing the **SET** button, you will move to the next setting.

#### 3. Select height


- Press the ▲ or ▼ button to select height (between 100 and 250 cm):

- Confirm your choice by pressing the **SET** button, you will move to the next setting.

#### 4. Select age

- Press the ▲ or ▼ button to select age (between 10 and 100 years):

#### 5. Save

- Press on the platform. The display shows "0.0 kg" and then the scale turns off after a few seconds.


### NOTES:


- When adjusting the settings, if no button is pressed for **6 seconds**, the parameters already pre-set will be saved automatically and the display will show "0.0 kg".
- If the display shows "0.0 kg" when adjusting the settings, press and hold the **SET** button for **3 seconds** to return to the setting mode and set the parameters (memory number, sex, size and age) again.

## Readings


- Press the **SET** button. The body fat scale will light up on the last selected memory.
- Use the **▲** and **▼** buttons to change the user (memory number).

- Then press your foot on the platform: when the display shows “0.0 kg”, you can step onto the platform. (**Note:** both your feet must be in contact with the electrodes otherwise the fat mass cannot be measured).

Then, the device takes the readings...  
... And displays the measured results 3 times and then saves them.


The body fat scale turns off after a few seconds.

The user's results will appear again as a reference for the next reading.

## Reminder of the results


- Press the **SET** button to turn on the body fat scale.

Wait 6 seconds or  
press on the platform


Press the button  
to switch to reminder mode

The screen successively displays the last saved user settings saved. Press **▲** to quickly move to the next setting.

The results scroll by 3 times and then the device turns off automatically.

## error messages

### Used battery indicator

Please replace the batteries with new ones with the exact same characteristics as the original ones.


### Overload indicator

The maximum weight accepted on the body fat scale has been exceeded. Please step off the platform to avoid damaging the unit.

### Indication of insufficient fat mass

The percentage of fat mass is too low.

### Indication of overly high fat mass

The percentage of fat mass is too high.

### Possible error

Please take the reading again to obtain a correct result.

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All information, designs, drawings and pictures in this document are the property of SOURCING & CREATION. SOURCING & CREATION reserves all rights to its brands, designs and information. Any copy and reproduction through any means shall be deemed and considered as counterfeiting.

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#### Protection of the environment

This symbol attached to the product means that it is an appliance whose disposal is subject to the directive on waste from electrical and electronic equipment (WEEE). This appliance may not in any way be treated as household waste and must be subject to a specific type of removal for this type of waste. Recycling and recovery systems are available in your area (waste removal) and by distributors. By taking your appliance at its end of life to a recycling facility, you will contribute to environmental conservation and prevent any harm to your health.

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#### \* Tested in our laboratories

Warranty valid from the date of purchase (receipt as proof of purchase). This warranty does not cover defects or damage caused by improper set up, incorrect use, or normal wear and tear of this product.

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