consignes d'usage

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instructions for use

CAREFULLY READ THESE INSTRUCTIONS BEFORE USING THE DEVICE AND KEEP THEM FOR FUTURE REFERENCE.

- This device is intended for household use only. It is not intended for any other use (professional, commercial, etc.) ...
- Only use the dental water flosser as indicated and for the purposes stated in the operating manual.
- Repairs cannot be performed by the user. Repairs must only be carried out by an authorised after-sales service.
- This device contains batteries that cannot be replaced.

- HydropulATORS may be used by children aged 8 and over and people with reduced physical, sensory or mental capabilities or insufficient experience or knowledge, provided they are supervised.
- Hydropumps should only be used with cold or hot water or specific solutions as defined in the manufacturer's instructions. Children must not play with the appliance.
- Children must be supervised if they are cleaning and maintaining the device.
- The charger should not be used outdoors or in damp places.
- Do not place the charger on or near a bathtub, shower, sink or any other container that holds water.
- Never submerge the charger in water or any other liquid.
- After cleaning the dental water flosser and before charging it, make sure that the charger is completely dry.
- Check that the voltage indicated on the charger's rating plate corresponds to the voltage of your electrical system.
- Only use the accessories provided with the device. Only use the charger provided with the device.
- Be sure to keep the device and charger away from hot surfaces.

- Never use this device if the power cord or charger is damaged, the device does not operate properly, you have dropped it, it is damaged, or it has fallen into water. Take it to an after-sales service centre for inspection and repair.
- For detailed information on recharging the battery, refer to the paragraph below in the manual.
- First connect the charger to the device, and then plug it into an easily accessible power point.
- Make sure that the power cord cannot be crushed or trapped, particularly near the power supply unit, the power point and where it connects to the device.
- Never unplug the device by pulling on the power cord.

- Always completely unroll the power cord when using it (to avoid overheating).
- Do not wrap the power cord around the charger when storing it.
- If the power cable is damaged, it must be replaced by the manufacturer, its after-sales service or similarly qualified persons in order to avoid danger.
- Always switch off the device (with the on/off button) when changing accessories and before cleaning.
- Unplug the charger after each use.
- Do not use on animals.
- For hygienic reasons, do not share your dental water flosser with other people.

- **WARNING:** if you have received oral or gum treatment, or if you have any oral problems, consult your dentist before using the device.
- If your gums bleed after using the device, consult your dentist.
- Do not exert excessive pressure when using the device.
- Use the device with cold or warm water only (max. 40°C). Risk of burns if hot water is used.
- Do not point the stream under your tongue or into your ears, eyes or other sensitive areas. Risk of injury!
- Do not clean the device with alcohol or other cleaning liquids.

Disposal of the battery

- The built-in rechargeable battery contains substances that can harm the environment.
- Before permanently disposing of this device, the rechargeable battery must be removed.
- The device must be disconnected from the mains supply when the battery is removed.
- The battery must be disposed of safely.
- Do not dispose of the battery or the device with household waste. To protect the environment, dispose of the device and batteries in accordance with current regulations.

- Place it in the collection bins provided for this purpose or take it back to your store.
- Never disassemble, crush or pierce the battery.
- Do not put it in a fire.
- Do not expose it to high temperatures or direct sunlight.
- Do not incinerate it.

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