## consignes d'usage

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## instructions for use

READ THESE INSTRUCTIONS CAREFULLY BEFORE USING THE APPLIANCE AND KEEP THEM FOR LATER REFERENCE AS NEEDED.

EN

- This appliance is intended for domestic use only. Any other use (commercial, professional, etc.) is not permitted.
- Use the water flosser only as indicated and for the uses described in the user guide.
- Use the electric toothbrush only for brushing your teeth. Any use other than that for which the appliance is intended must

be considered improper and therefore dangerous.

- The user may not perform any repairs. All repairs must be carried out exclusively by an authorised after-sales service department.
- The toothbrushes can be used by children and adults whose physical, sensory or mental capabilities are impaired or whose experience or knowledge is insufficient, provided they are supervised or instructed on how to use them safely, and provided they understand the potential dangers of doing so. The device must not be cleaned or maintained by children without supervision.

- The water flossers can be used by children over the age of 8 and adults whose physical, sensory or mental capabilities are impaired or whose experience or knowledge is insufficient, provided they are supervised.
- Children must not play with this device.
- The charging base should not be used outdoors or in damp places.
- Do not place the charging base on or near a bath, shower, washbasin or near any other container filled with water.
- Never immerse the charging base in water or any other liquid.

- After cleaning the water flosser or toothbrush and before charging these devices, make sure the charging base is completely dry.
- Check that the supply voltage on the charging base data plate is compatible with your electrical installation.
- Only use the accessories supplied with this device. Use only the charging base supplied with the device.
- Keep the device and charging base away from hot surfaces.
- Never use this device if the power cord or charging base are damaged, if the device does not work properly, if it has been

dropped, if it is damaged or has been immersed in water. Please take it to a customer service department centre, where it will be inspected and repaired.

- The battery cannot be replaced by the user, please go to your Customer Service Centre if it needs to be replaced.
- First plug the mains charger into the device and then plug the mains charger into an easily accessible power outlet.
- Make sure the power cord cannot be crushed or trapped, especially around the adapter, power outlet and where it is connected to the appliance.

- Never move the device by pulling on the power cord.
- Always unwind the power cord completely when in use (to prevent it from overheating).
- Do not wrap the power cord around the charging base when storing it.
- Use only the charging base supplied with the device. Do not use the charging base for other purposes.
- The power cable cannot be replaced. If the cable is damaged, the power supply should be discarded. Contact your store to obtain replacement parts.
- Always switch off the device (using the On/Off button) when

changing accessories (brush heads or nozzles) and before cleaning.

- Unplug the charging station after each use.
- Do not use on animals.
- For hygiene reasons, do not lend your water flosser or toothbrush to anyone else.
- **WARNING:** if you have had oral work done or treatment for gingivitis, or if you suffer from oral problems, ask your dentist for advice before using the device.
- If your gums bleed after using the device, consult your dentist.
- Do not apply excessive force when using the device.

- Use the device only with cold or warm water (40°C max.). Risk of burns if hot water is used.
- Do not direct the spray under the tongue or in the ears, eyes or any other sensitive areas. Risk of injury!
- Do not clean the device with alcohol or any other liquid detergents.

### **Battery waste separation**

- The built-in rechargeable batteries contain substances which can harm the environment.
- Do not dispose of the battery or appliance with household waste.

- In order to protect the environment, the battery and the appliance must be disposed of in accordance with applicable regulations.
- Dispose of the appliance in the appropriate recycling container or return it to your retail store.
- Never disassemble, crush or pierce the battery.
- Do not dispose of the battery in a fire. Do not burn it.
- Do not expose it to high temperatures or direct sunlight.

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