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How to use

Before using your Bluetens is turned off before removing electrodes.

Electrodes do not stick well if any lotion, oil, etc is left on the skin. Clean and dry the skin prior to each application.

After the cable with the electrodes prior to application. The electrodes following the positioning proposed on the App. Press firmly on the centre of the electrodes to the outer surface of the electrodes must be in full contact with your skin. Adhesion improves when the electrodes reach skin temperature.

After the electrodes back on the right side of the protective plastic liner after use and put the liner in the sealed bag to prevent gel drying. Take out the electrodes by pulling on the corner of the electrode and slowly lift it from the skin, touching the adhesive gel as you lift it.

little as possible.

6) The electrodes are intended for single patient use only. The number of application depends on your skin condition, anatomical placement, type of treatment and care after use. If adhesive gel becomes dry you can add one or two drops of ordinary tap water spread evenly and gently with your fingertips. Do not use too much water as this may reduce the life of the gel.

7) Store the electrodes bag at room or cool temperature (+5 - +25 °C) and keep out from direct sunlight.

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