

| TÊTE DE BROSSE | | |
|---|-----------------------------|--|
|  | Nom du produit | Tête de brosse compatible MyVariations® |
| | Product name | Brush head |
| | EAN | S : 3770031226041 M : 3770031226034 |
| | SKU | S : Y2BshBrushUnivNaWtWtS M : Y2BshBrushUnivNaWtWtM |
| | Catégorie / Category | Têtes de brosse / Brush heads |

| INFORMATIONS PRODUIT FR | | | |
|---|---|--------------------|---------------------------------|
| Description courte | Description longue | Contenu | Visuels |
| <p>Les têtes de brosse Y-Brush, compatibles avec Philips® Sonicare®, sont flexibles pour s'adapter à toutes les mâchoires et épouser l'arcade dentaire. Entièrement tapissées de poils en nylon médical, elles permettent de brosser toutes les faces de toutes les dents !</p> | <p>Quand changer sa tête de brosse ? Tous les 2 à 6 mois, dès que vous voyez que les poils sont abîmés.</p> <p>Les têtes de brosse Y-Brush sont conçues sans latex ni silicone. Elles disposent de 35 000 filaments en nylon médical souples afin de retirer efficacement la plaque dentaire. Ces poils sont inclinés à 45° afin de respecter la méthode de Bass.</p> <p>Guide des tailles et produits compatibles Taille adulte (M) Pour toute personne âgée de 12 ans ou plus.</p> <p>Compatible avec les différents modèles de brosses à dents de la marque MyVariations®.</p> | - 1 tête de brosse | Têtes de brosse |

| INFORMATIONS PRODUIT EN | | | |
|---|--|----------------|---------------------------------|
| Short description | Long description | Content | Visuals |
| <p>Y-Brush brush heads, compatibles with Philips® Sonicare®, are flexible to adapt to all jaw sizes and fit the dental arch. Completely covered with medical-grade nylon bristles, they allow brushing all surfaces of all teeth!</p> | <p>When to change the brush head? Every 2 to 6 months, as soon as you see that the bristles are worn out. Y-Brush brush heads are designed without latex or silicone. They have 35,000 soft medical-grade nylon filaments to effectively remove dental plaque. These bristles are angled at 45° to follow the Bass method.</p> <p>Size guide and Compatible products Adult Size (M) For everyone aged 12 and above.</p> <p>Compatible with MyVariations® Toothbrush Models.</p> <p>This specific product DOES NOT FIT Y-Brush handles (NylonStart, NylonKids, NylonBlack).</p> <p>*Important Disclaimer: MyVariations® is registered trademark of Le</p> | - 1 brush head | Têtes de brosse |

| INFORMATION MARKETING / MARKETING INFORMATION | | | |
|---|---|---|--|
| Claim principal / positionnement produit | Passez facilement du brossage sonique traditionnel au brossage Y-Brush à l'aide d'un seul manche MyVariations®. Deux têtes, un manche, une commodité illimitée. | Main claim / product positioning | Switch easily from traditional sonic brushing to Y-Brush brushing with just one MyVariations® handle. Two heads, one handle, unlimited |
| Spécificités techniques | <ul style="list-style-type: none"> - Super flexible - Sans latex ni silicone - 35 000 filaments en nylon souple - Inclinaison des poils à 45° - Testée cliniquement | Technical features | <ul style="list-style-type: none"> - Super flexible - Latex- and silicone-free - 35,000 soft nylon filaments - 45° hair inclination - Clinically tested |
| Bénéfices conso | Confortable, toutes les faces de toutes les dents brossées, méthode de brossage recommandée suivie | Consumer benefits | Comfortable, all sides of all teeth brushed, recommended brushing method followed |
| Entretien | <p>Rinçage de la brosse après chaque utilisation, passage du doigt entre les poils si nécessaire.</p> <p>Possibilité également de la faire tremper dans du bain de bouche toute une nuit, 1 à 2 fois par semaine : 1/3 de bain de bouche + 2/3 d'eau.</p> | Maintenance | <p>Rinse brush after each use, running your finger between the bristles if necessary.</p> <p>You can also soak the brush in mouthwash overnight, once or twice a week: 1/3 mouthwash + 2/3 water.</p> |