# beurer

Phototherapy unit

TL 90

<b>D</b>	<b>Tageslichtlampe</b>	
	Gebrauchsanweisung .....	2–8
<b>GB</b>	<b>Daylight therapy lamp</b>	
	Instruction for use .....	9–14
<b>F</b>	<b>Lampe de luminothérapie</b>	
	Mode d'emploi .....	15–20
<b>S</b>	<b>Dagsljuslampa</b>	
	Bruksanvisning .....	21–26
<b>FIN</b>	<b>Kirkasvalolamppu</b>	
	Käyttöohje .....	27–32
<b>DK</b>	<b>Dagslyslampe</b>	
	Brugsanvisning .....	33–38

<b>N</b>	<b>Dagslyslampe</b>	
	Bruksanvisning .....	39–44
<b>TR</b>	<b>Gün ışığı lambası</b>	
	Kullanma Talimatı .....	45–50
<b>RUS</b>	<b>Прибор дневного света</b>	
	Инструкция по применению .....	51–57
<b>PL</b>	<b>Lampa o świetle dziennym</b>	
	Instrukcja obsługi .....	58–63

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Germany

## Content

1. About the daylight therapy lamp ..... 9	7. Cleaning and care of the unit ..... 13
2. Signs and symbols ..... 9	8. Replacing the lamp tubes /starter..... 13
3. Notes ..... 10	9. Storage ..... 13
4. Unit description ..... 11	10. Disposal..... 13
5. Initial use ..... 11	11. What if there are problems? ..... 13
6. Operation..... 11	12. Technical specifications ..... 14

## 1. About the daylight therapy lamp

### Dear valued customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air.

With kind regards  
Your Beurer team

### Why use a daylight therapy lamp?

When the hours of daylight are noticeably decreased in autumn and people increasingly stay inside, the effects of a lack of light may cumulate. This is often described as „winter depression“. The symptoms can present themselves in a number of ways:

- Imbalance
- Subdued mood
- Lack of energy and listlessness
- Generally feeling under the weather
- Need for more sleep
- Loss of appetite
- Difficulty concentrating

The cause of these symptoms is the fact that light, particularly sunlight, is essential for life and has a direct effect upon the human body. Sunlight indirectly controls the production of melatonin, which is released into the blood only in darkness. This hormone tells the body when it is time to sleep. In the months when there is less sunlight, therefore, production of melatonin increases. This makes it harder to get up in the morning, since the body's functions are shut down. If you use the daylight therapy lamp immediately after waking up in the morning—thus as early as possible—you can end the production of melatonin so that your mood takes a positive turn. In addition, light deficiency also hinders production of the happiness hormone serotonin, which decisively affects our well-being. So application of light involves quantitative changes in hormones and neurotransmitters in the brain, which can affect our level of activity, our feelings and our well-being. Light can help reverse this hormone imbalance. To compensate for such a hormonal imbalance, the daylight therapy lamp can create a suitable replacement for natural sunlight.

In the medical field, daylight therapy lamps are used to compensate for the effects of a lack of light. Daylight therapy lamps simulate daylight over 10,000 lux. The light can influence the human body and be used as a treatment or as a preventative measure. Normal electric light is not sufficient to influence the hormonal balance. Because in a well-lit office the illuminance level amounts to just 500 lux.

- Included in delivery:**
- Daylight therapy lamp
  - These operating instructions

## 2. Signs and symbols

The following symbols appear in these instructions or on the device.

	Read the instructions of use	 <b>Note</b>	Note on important information
 <b>Warning</b>	Warning instruction indicating a risk of injury or damage to health		Manufacturer

 <b>Important</b>	Safety note indicating possible damage to the unit/accessory		Device in protection class 2
	Serial number		OFF / ON
	The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products	<b>IP21</b>	Protected against solid foreign objects 12.5 mm in diameter and larger, and against vertically falling drops of water
<b>Storage</b> 	Permissible storage temperature and humidity	<b>Operating</b> 	Permissible operating temperature and humidity

### 3. Notes

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

#### Safety notes

##### Warning

- The light of the daylight therapy lamp is only intended for use on the human body.
- Before use, ensure that all packaging materials are removed and that there is no visible damage to the unit or accessories. When in doubt, do not use the unit and contact your dealer or the customer service address provided.
- Also ensure that the daylight therapy lamp has a firm footing.
- The mandatory minimum distance between the daylight therapy lamp and flammable objects is 1.0 m.
- Connect the unit only to the mains voltage listed on the type plate.
- Do not dip the unit into water and do not use it in wet rooms.
- Do not use it in the vicinity of people who are helpless, children under age 3 or people who are insensitive to heat (people with skin lesions due to illness); do not use after taking painkillers or drugs or drinking alcohol.
- Keep children away from packaging materials (risk of suffocation).
- Do not cover up or pack away the unit while it is warm.
- Always remove the mains plug and let the unit cool off before you touch it.
- Do not touch the unit with wet hands while it is plugged in; do not allow any water to be sprayed onto the unit. The unit must be operated only when it is completely dry.
- Keep the mains cable away from hot objects and naked flames.
- Protect the unit from strong impacts.
- Do not pull the mains plug from the socket by the cord.
- If the mains connection cable of this device is damaged, it must be disposed of. If it cannot be removed, the device must be disposed of.
- Disconnection from the mains supply is only guaranteed if the mains plug has been removed from the socket.
- Do not use the device in the presence of flammable anaesthetic gas connections with air, oxygen or nitrogen oxide.
- No calibration and no preventative checks or maintenance is to be carried out on this device.
- You cannot repair the device. The device contains no parts that you can repair.
- Do not make any changes to the device without the manufacturer's permission.
- If the device has been changed, thorough tests and checks must be carried out to ensure the continued safety of further use of the device.
- As the patient/user, you can operate all functions of this device.

## General notes

### ⚠ Important

- If you are taking medication such as pain relievers, antihypertensives or antidepressants, consult your doctor before using the daylight therapy lamp.
- People with retinal diseases as well as diabetics need to be examined by an ophthalmologist before beginning to use the daylight therapy lamp.
- Please do not use if you have an eye disease such as cataracts, glaucoma, diseases of the optic nerve in general and inflammation of the vitreous humour.
- If you have health concerns of any kind, contact your general practitioner!
- Before you use the unit for the first time, remove all packaging materials.
- The lamp tubes are excluded from the warranty.

## Repairs


### ⚠ Important

- Do not open the device. Please do not attempt to repair the unit yourself. This could result in serious injuries. Failure to observe this regulation shall void the warranty.
- If you need to have the unit repaired, contact customer service or an authorised dealer.

## 4. Unit description

### Overview

1. Fluorescent screen
2. Rear of the device housing
3. Stand
4. On/Off button
5. LED display

## 5. Initial use

Take the device out of the plastic wrapping. Check the device for damage or faults. If you notice any damage or faults on the device, do not use it and contact customer services or your supplier.

### Positioning the device

Place the device on an even surface. The position should be chosen to ensure a distance of between 15 cm and 50 cm between the user and the device. The lamp is most effective at this distance.

### Mains connection

Connect the unit only to the mains voltage listed on the type plate. Push the plug completely into the mains socket.

### ⚠ Note:

- Ensure that there is a socket near the place where the lamp will stand.
- Lay the mains cable in such a way that no one can trip over it.

## 6. Operation

### 1 Switching on the lamp

▶ Press the On/Off button. The tubes turn on. Set the angle of the fluorescent screen in accordance with your individual requirements.

### ATTENTION!

The brightness of the tubes may be visibly different when switched on. The tubes develop their full luminosity after a couple of minutes.

2	<p><b>Enjoying the light</b></p> <p>Sit as close as possible to the lamp, between 15 cm and 50 cm. You can go about your normal activities while using it. You can read, write, make telephone calls, etc.</p> <ul style="list-style-type: none"> <li>• Every so often, briefly look directly into the light, since the effect results from the eyes/ retina receiving the light.</li> <li>• Use the daylight therapy lamp as often as you want. However, the treatment is at its most effective if you carry out the light therapy according to the prescribed time period of at least 7 successive days.</li> <li>• The most effective time of day for the treatment is between 6am and 8pm and we recommend that it is used for 2 hours per day.</li> <li>• Do not look directly into the light the whole time you are using it, however, since that could lead to overstimulation of the retina.</li> <li>• Start with shorter periods of illumination, which you increase over the course of a week.</li> </ul> <p><b>Note:</b></p> <p>Eye aches and headaches may arise after the first times using the lamp. These will go away in later sessions, as the nervous system will become accustomed to the new stimulation.</p>																				
3	<p><b>LED display</b></p> <p>The daylight therapy lamp TL90 indicates your current treatment time with the aid of 8 LEDs. When the device is switched on, all LEDs will briefly light up. After a treatment time of 15 minutes, one LED will light up at a time.</p> <table border="1" data-bbox="141 600 787 914"> <thead> <tr> <th>Treatment time</th> <th>Number of lit LEDs</th> </tr> </thead> <tbody> <tr> <td>0 – 15 minutes</td> <td>0</td> </tr> <tr> <td>16 – 30 minutes</td> <td>1</td> </tr> <tr> <td>31 – 45 minutes</td> <td>2</td> </tr> <tr> <td>46 – 60 minutes</td> <td>3</td> </tr> <tr> <td>61 – 75 minutes</td> <td>4</td> </tr> <tr> <td>76 – 90 minutes</td> <td>5</td> </tr> <tr> <td>91 – 105 minutes</td> <td>6</td> </tr> <tr> <td>106 – 120 minutes</td> <td>7</td> </tr> <tr> <td>&gt;120 minutes</td> <td>8</td> </tr> </tbody> </table> <p>The maximum recommended treatment time will pass after 120 minutes and all LEDs will begin to flash 5 times. After the 120 minutes has passed, switch off the daylight therapy lamp by pressing the ON/OFF button.</p>	Treatment time	Number of lit LEDs	0 – 15 minutes	0	16 – 30 minutes	1	31 – 45 minutes	2	46 – 60 minutes	3	61 – 75 minutes	4	76 – 90 minutes	5	91 – 105 minutes	6	106 – 120 minutes	7	>120 minutes	8
Treatment time	Number of lit LEDs																				
0 – 15 minutes	0																				
16 – 30 minutes	1																				
31 – 45 minutes	2																				
46 – 60 minutes	3																				
61 – 75 minutes	4																				
76 – 90 minutes	5																				
91 – 105 minutes	6																				
106 – 120 minutes	7																				
>120 minutes	8																				
4	<p><b>Important instructions</b></p> <p>When using the lamp, maintain the recommended distance of 15 – 50 cm between your face and the lamp. The duration of the application depends on the distance:</p> <table border="1" data-bbox="141 1121 868 1268"> <thead> <tr> <th>Lux</th> <th>Distance</th> <th>Application duration</th> </tr> </thead> <tbody> <tr> <td>10.000</td> <td>approx. 15 cm</td> <td>0,5 hour</td> </tr> <tr> <td>5.000</td> <td>approx. 30 cm</td> <td>1 hour</td> </tr> <tr> <td>2.500</td> <td>approx. 50 cm</td> <td>2 hours</td> </tr> </tbody> </table> <p>Basically, the closer you are to the source of light, the less time you should use it.</p>	Lux	Distance	Application duration	10.000	approx. 15 cm	0,5 hour	5.000	approx. 30 cm	1 hour	2.500	approx. 50 cm	2 hours								
Lux	Distance	Application duration																			
10.000	approx. 15 cm	0,5 hour																			
5.000	approx. 30 cm	1 hour																			
2.500	approx. 50 cm	2 hours																			
5	<p><b>Enjoying light over longer periods</b></p> <p>Repeat your use of the lamp in the dark seasons for at least 7 successive days, or even longer, depending on your individual needs. If possible, conduct the treatment in the morning hours.</p>																				

<b>6</b>	<p><b>Switching off the lamp</b></p> <p>▶ Press the On/Off button. The tubes turn off. Pull the mains plug from the socket.</p> <p><b>CAUTION!</b></p> <p>The lamp is hot after being used. Let the lamp cool off first long enough before you put it away and/or pack it up!</p>
----------	---

## 7. Cleaning and care of the unit

The unit should be cleaned occasionally.

### Important

- Ensure that no water penetrates inside the unit!  
Before you clean the unit, always switch it off, unplug it and let it cool off.
- Do not wash the unit in a washing machine!  
Clean using a damp cloth, to which you can apply a little detergent if necessary.
- Do not touch the unit with wet hands while it is plugged in; do not allow any water to be sprayed onto the unit. The unit must be operated only when it is completely dry.  
Do not use abrasive detergents and never immerse the unit in water.

## 8. Replacing the lamp tubes/starter

The service life of the fluorescent tubes is on average 10,000 hours. If one tube needs replacing, both fluorescent tubes must be replaced with new tubes of the same type. For this purpose, contact a retailer, an authorised service centre or appropriately qualified personnel or contact the customer service.

The following tubes are approved as replacement:

- Philips – PL-L36W/865/4P
- The manufacturer’s tubes are available with the order number 162.942 from customer support

### Warning

- For environmental reasons, do not dispose of the fluorescent tubes with the household waste. Dispose of the unit at a suitable local collection or recycling point. Please contact the local authorities responsible for waste disposal.
- Unplug the device, to avoid the risk of electrocution.
- Take care when handling fluorescent tubes to avoid the risk of cut injuries.

## 9. Storage

If you are not going to use the device for an extended period of time, disconnect it and store it in a dry place, out of the reach of children.

Follow the storage instructions provided in the “Technical data” chapter.

## 10. Disposal

For environmental reasons, do not dispose of the device in the household waste at the end of its useful life. Dispose of the device at a suitable local collection or recycling point. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

## 11. What if there are problems?

Problem	Possible cause	Solution
Device does not light up	On/off button switched off	Switch on the on/off button.
	No electricity	Connect the mains plug correctly.
	No electricity	The mains cable is defective. Contact customer services or your retailer.
	Fluorescent tube has reached the end of its service life. Fluorescent tube is defective.	If you need to have the unit repaired, contact customer service or an authorised dealer.

## 12. Technical specifications

Dimensions (WHD)	337 x 581 x 241 mm
Weight	~ 3 kg
Lighting elements	FSL-YDW-36HRR
Power	2 x 36 Watts
Light intensity	10.000 lux (Application approx. 15 cm)
Radiation	Output of radiation beyond the visible spectrum (infrared and UV) is so low that it is harmless to eyes and skin.
Mains connection	230 V ~, 50 Hz
Operating conditions	+5°C to +45°C, 30 – 90 % relative humidity; 700 – 1060 hPa ambient pressure
Storage conditions	+5°C to +50°C, 30 – 90 % relative humidity; 700 – 1060 hPa ambient pressure
Product classification	Protection class II
Degree of efficiency of the ballast	~ 98%
Source efficiency (lm/W) of the bulb	64 lm/W
Color temperature of tubes	6.500 Kelvin
Fuse type	F2AL, 250V

Subject to technical modifications.

The serial number is located on the device or in the battery compartment.

This device meets the requirements of the European Directive 93/42/EEC and 2007/47/EEC Medical Devices as well as the Medical Devices Act.

This device complies with European Standard EN60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility. Please note that portable and mobile HF communication systems may interfere with this unit. More details can be requested from the stated Customer Service address or found at the end of the instructions for use.

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# ELECTROMAGNETIC COMPATIBILITY INFORMATION

<b>Guidance &amp; Declaration – electromagnetic immunity</b>			
The model TL90 is intended for use in the electromagnetic environment specified below. The customer or the user of the model TL90 should assure that It is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment – guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0,5 cycle  40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles  70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles  <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 sec	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0,5 cycle  40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles  70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles  <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the TL90 requires continued operation during power mains interruptions, it is recommended that the TL90 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<b>NOTE</b> $U_T$ is the a.c. mains voltage prior to application of the test level.			

**Guidance & Declaration – electromagnetic immunity**

The TL90 is intended for use in the electromagnetic environment specified below. The customer or the user of the TL90 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the TL90, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> <p>3V</p> <p><math>d=1.2 P^{1/2}</math> 80 MHz to 800 MHz</p> <p><math>d=2.3 P^{1/2}</math> 800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following</p> 

**NOTE 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the TL90 is used exceeds the applicable RF compliance level above, the daylight therapy lamp should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the TL90.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

**Recommended separation distances between  
portable and mobile RF communications equipment and the TL90**

The TL90 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the daylight therapy lamp can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the TL90 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \times P^{1/2}$	80 MHz to 800 MHz $d = 1.2 \times P^{1/2}$	800 MHz to 2,5 GHz $d = 2.3 \times P^{1/2}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Taishan Anson Electrical Appliances Co. Ltd • Dragon Mountain Industrial Estate  
Duanfen • Taishan, Guangdong Province, China

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Distributed by: BEURER GmbH • Söflinger Str. 218 • 89077 Ulm, Germany
designed4inspiration GmbH • Feldgasse 34 • A-7072 Mörbisch am See