 More detailed product information is available at www.slendertone.com. This includes ab belt setup information, how to position the belt properly for optimal results and how to use the range of programmes available to best suit your individual needs.

 Use this product only for its intended purpose. For best results and maximum comfort, we recommend that you use genuine Slendertone accessories, do not use gel pads or other materials not recommended by Slendertone.

PRODUCT WARRANTY

Should your Slendertone product develop a fault within two years of purchase, Slendertone will undertake to replace or repair the product or any parts found to be defective with no charge for labour or materials*, provided the product:

- Has been used for its intended purpose.
- Has not been connected to an unsuitable power source.
- Has not been subjected to misuse or neglect.
- Has not been modified or repaired by anyone other than an approved Slendertone agent.

This guarantee complements existing national guarantee obligations and does not affect your statutory rights as a consumer.

* This excludes consumables (e.g. gel pads, belt, etc.) which are subject to normal wear and tear.

CONTRA-INDICATIONS, PRECAUTIONS & WARNINGS

 Your Slendertone product is suitable for use by all healthy adults. However, as with other forms of exercise, some care is needed when using it, so always follow the points below and read your user manual carefully before use. Some of the points below are gender specific.

Please do not use if:

- You have an electronic implant (e.g. cardiac pacemaker or defibrillator) or suffer from any other heart problem.
- You are pregnant.
- You suffer from cancer, epilepsy or are under medical supervision for cognitive dysfunction.
- The controller is in close proximity (e.g. 1m) to shortwave or microwave therapy equipment.
- You are connected to high-frequency surgical equipment.
- Wearing the device necessitates placement over areas at which drugs/medicines are administered by injection (short term or long term) e.g. hormone treatment.

Please wait before using your Slendertone product until:

- At least six weeks after the birth of your baby (you must consult your doctor first).
- At least three months after having a Caesarean (you must consult your doctor first).
- One month after an IUD contraceptive device (e.g. coil) has been fitted.
- The heavy days of your period have finished.

Please get your doctor's advice before using your Slendertone product if:

- You wish to use a Slendertone belt and have a bad back. You should also ensure the intensity is kept low.
- You have any serious illness or injury not mentioned in this guide.
- You have recently had an operation.
- You take insulin for diabetes.
- You want to use it on a young child.
- You suffer from muscle or joint problems.

When applying the gel pads and belt, always remember to:

- Place the gel pads and belt ONLY on the abdomen. For more information, please see instructional video on Slendertone.com/evolve.
- Never place the gel pads on the front or sides of the neck, across or through the heart (i.e. one gel pad on the front of the chest and one on the back), in the genital region or on the head.
- Application of gel pads near the thorax may increase the risk of cardiac fibrillation.
- Avoid any recent scars, broken or inflamed skin, areas of infection or susceptibility to acne, thrombosis or other vascular problems (e.g. varicose veins), or any parts of the body where feeling is limited.
- Avoid areas of injury or restricted movement (e.g. fractures or sprains).

The controller contains 1 x 3.7 V Nominal Lithium Ion Polymer Battery.

- The Output Frequency indicates the number of pulses per second transmitted by the controller.
- This is measured in hertz, which is indicated by "Hz".
- Output (RMSA) means the maximum output root mean square current for each channel.
- Output (RMSV) means the maximum output root mean square voltage for each channel.

Attention

This product can deliver currents in excess of 2mA/cm² if used with incorrect electrodes.

Rated Outputs Voltage/Currents: Type 735

Parameter	500Ω	1 KΩ	1.5 kΩ
Output RMSV	8.4V	14.4V	15.5V
Output RMSA	16.7mA	14.4mA	10.3mA
Output Frequency	50-75Hz	50-75Hz	50-75Hz
DC Component: Approx.	0 C	0 C	0 C
Positive Pulse Width	200–350µs	200–350µs	200–350µs
Negative Pulse Width	200–350µs	200–350µs	200–350µs
Interphase Interval	100µs	100µs	100µs

- Avoid placing the gel pads directly over metal implants.
- Do not switch the controller on until it has been correctly inserted into the holster on the belt.

Possible adverse reactions:

- A small number of isolated skin reactions have been reported by people using muscle stimulation devices, including allergies, a prolonged reddening of the skin and acne.

To reposition pads during a session:

- Always pause the programme currently running, unfasten the belt and then refasten it behind your back once the position of the pads has been adjusted.
- After strenuous exercise or exertion always use a lower intensity to avoid muscle fatigue.

General safekeeping:

- Keep your product out of the reach of children.
- The studs and gel pads must not be connected to other objects.
- Do not use your controller at the same time as any other device which transfers an electrical current into the body (e.g. another muscle stimulator).
- Cease using your product if you are feeling lightheaded or faint. Consult your doctor if this happens.
- Do not touch the gel pads or metal studs while the controller is switched on.
- Do not use while driving, operating machinery or cycling.
- For first time users, muscle stimulation can be an unusual sensation. We recommend that you begin with the initiation programme at low stimulation intensities to familiarise yourself with the sensation before progressing to higher stimulation intensities.
- In rare situations, first time users may feel faint. We recommend that you use the product while seated or lying down until you become accustomed to the sensation.
- For hygiene reasons the belt is for use only by one person.
- Do not share your belt with anyone else.
- Do not overexert yourself when using muscle stimulation. Any workout should be at a level comfortable for you.
- Do not use your product if you are wearing a belly-button ring. Remove the ring before you begin a session.

- Slendertone will not accept responsibility if the guidelines and instructions supplied with this product are not followed.
- Although compliant with applicable EMC requirements, this device may still interfere with more sensitive equipment, please move away or switch off.
- Keep this device out of the reach of pets or pests.
- Always allow two hours for the controller to cool down between stimulation sessions. During normal use, the controller generates a small amount of heat and tests have shown that in extreme conditions this can cause a temperature rise up to 3°C / 5.4°F under the belt at the controller location. At high ambient temperatures, close to 40°C / 104°F, care is needed to avoid prolonged exposure to this increased temperature.

-  If you are in any doubt about using your Slendertone belt for any reason, please consult your doctor before use.

 For a full list of frequently asked questions, visit www.slendertone.com

BATTERY

Battery Charging & Discharging

- Your device is powered by a rechargeable battery. The battery can be charged and discharged hundreds of times, but it will eventually wear out. To ensure optimum battery performance you should use your device regularly.
- To conserve energy, disconnect the USB cable from the controller when the battery is fully charged.
- Leaving the battery in hot or cold places, such as in a closed car in summer or winter conditions, may reduce the capacity and lifetime of the battery. A device with a hot or cold battery may cease working temporarily, even when the battery is full charged. Battery performance is particularly limited in temperatures below freezing.
- Do not dispose of batteries in a fire as they may explode. Batteries may also explode if damaged. Dispose of batteries according to local regulations. Please recycle when possible. Do not dispose as household waste.

Charging Specifications

When charging the battery, you must use a charger with the following outputs:

Output Voltage: 5V DC nominal.
Output Current: 0.5A/500mA minimum.
Output Connector: Micro-B USB Connector.

PRODUCT CARE

Caring For Your Belt

Your belt may be cleaned using a lightly dampened cloth, but you must first remove the controller and gel pads. Always follow the instructions on the label when washing the belt.

-  Never machine-wash or hand-wash your belt. Always clean it using a lightly dampened cloth or sponge.

-  Do not use bleach when cleaning the belt.

-  Do not dry clean your belt.

-  Do not tumble dry your belt. Dry the belt on a flat surface. Do not dry it over anything hot (e.g. a radiator) as the belt contains plastic parts. Ensure the belt is completely dry before using it again.

-  The belt should not be ironed.

Belt Materials:
Outer material: 87% Nylon, 13% Elastane;
Hook & Loop: 100% Nylon;
Inner Foam: 100% Polyurethane;
Textured Panel: 70% Polyurethane, 30% Polyester.

Caring For Your Controller

Your controller should not be allowed to get wet or be left in excessive sunlight. It may be cleaned using a soft, lightly dampened cloth. Do not allow the interior of your controller to become wet. Do not use detergents, alcohol, spray aerosols or strong solvents on your controller. If your controller is damaged, you should not use it, but you should contact the Slendertone Customer Care Helpline for advice. Repairs, service and modifications may not be carried out by anyone other than qualified service personnel authorised by Slendertone.

-  Access to the interior of the controller is not required for maintenance purposes. Do not open your controller casing or alter the product in any way that could cause harm to individuals. Tampering with the product would also invalidate the warranty.

Caring For Your Gel Pads

The gel pads are for a single person use only. You should ensure your skin is clean and free from oils, creams and other lotions before use. Avoid touching the black side of the gel pads with your fingers as much as possible. The durability and effectiveness of the gel pads depends entirely on proper usage, storage and care.

After use, replace the covers on the black side of the gel pads and store your product in a cool, dry place until the next use.

Gel Pads should be replaced after 20-30 toning sessions.

-  Do not use gel pads that are damaged in any way.

-  For a full list of frequently asked questions, visit www.slendertone.com

ACCESSORIES

Under no circumstances should anything other than Slendertone accessories be used with your Slendertone product (type 735). Any others may not be compatible with your product and could degrade the minimum safety levels. You can purchase all accessories from www.slendertone.com or by calling customer care.

- 1 large adhesive gel pad (Type 732)
- 2 small adhesive gel pads (Type 733)

USB Cable: Type A to Micro B USB Charger

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