
Programme Guide

SLENDERTONE®

Evolve Abs

Guide des Programmes
Guía del Programa
Guida ai Programmi

- Programmübersicht
- Programmahandleiding
- Guia do Programa

ENG Welcome

Congratulations on purchasing your Slendertone Evolve Abs! You are a part of a growing, supportive community of #Slendertoners who want to look and feel better and who use our ab toning belt to help them achieve this goal!

With 12 Toning Programmes to select from, including our new 6-week progressive toning plan and 15-minute "Abs Blitz" feature, the Evolve Abs has something to suit everybody. So, let's get started!

To Set-Up your Evolve Abs, Scan the QR Code above or type the URL into your web browser: <http://www.slendertone.com/pages/evolve-abs-set-up>

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Programmes

With 12 passive or active toning programmes to choose from, including a 6-week progressive toning plan, our powerful Abs Optimiser programmes and new Abs Blitz feature, there is something to suit everybody!

PASSIVE TONING

(i.e. only the ab belt is working to stimulate the muscles)

Classic Toning

Prog.	Description	Duration
1	Initiation Perfect for first-time Electrical Muscle Stimulation (EMS) users. Familiarise yourself with the unique EMS sensation before moving to the more challenging intermediate level.	20 Mins
2	Intermediate For those who want to build on their beginners' level of experience with EMS, or those coming back to EMS training after an extended break. Try this programme before embarking on the 6 week training plan.	30 Mins

Progressive Toning Plan

Prog.	Description	Duration
3	Weeks 1-2 Weeks 1-2 are the "base camp" of the Progressive Toning Plan. Complete 5 toning sessions each week and proceed to the more challenging plan in weeks 3-4.	30 Mins
4	Weeks 3-4 Keep climbing! You are now ready to progress to the more challenging weeks 3-4. Don't forget to keep increasing the intensity. Complete 5 toning sessions each week and you will begin to see real results.	30 Mins
5	Weeks 5-6 Almost there! Weeks 5-6 are an advanced EMS workout delivered over 30 minutes. Increase the intensity to a comfortable but challenging level and complete 5 times a week to reach the summit of your toning plan!	30 Mins

Abs Optimiser

Prog.	Description	Duration
6	Ab Power The most powerful programme - characterised by multiple cycles of rest and workout phases - delivers a high repetition total Abs workout. Use this once you have completed your 6-week programme to maintain and build on your newly acquired abdominal strength and definition!	30 Mins
7	Ab Power Plus Looking to really challenge those abs for ultimate definition? This Advanced Abs Activator, the longest toning session at 45 minutes, will give you a powerful, strong toning session.	45 Mins

Abs Blitz

Prog.	Description	Duration
8	Total Abs Blitz A top-up or stand-alone workout, the Abs Blitz is a challenging and targeted 15 minute intense session which delivers multiple phases of workout activation for more toned and firmer abs.	15 Mins
9	Oblique Abs Blitz A top-up or stand-alone workout, this 15 minute session focuses on the oblique or "love handle" ab muscles. Delivering multiple pulse parameters, this workout focuses on firming and strengthening these shape-enhancing muscles.	15 Mins
10	Power up your Plank Use this programme to build the strength required to carry out one of the most effective ab muscle exercises around - the Plank. This 10 minute session is specifically designed to build muscle strength and endurance to assist and increase your plank hold, but is not designed for use while completing the Plank exercise.	10 Mins

ACTIVE TONING

(i.e. wear the ab belt while carrying out the following light exercises to stimulate your ab muscles)

Prog.	Description	Duration
11	Cardio Toning Maximise the benefits of any light cardio toning activity such as walking, light jogging or HIIT workouts, to strengthen the abdominal muscles while you work on your cardio resistance.	40 Mins
12	Plank Partner Designed to be undertaken while doing the Plank, these short bursts are specifically designed to complement the plank exercise and develop muscle strength and endurance to assist your plank hold.	2 Min Bursts

Note:

The success of your toning sessions depends on the quality of your Gel Pads, which should be replaced after 25-30 sessions. Gel Pads can be purchased on slendertone.com

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Copyright © 2020 Bio-Medical Research Ltd. All rights reserved.
SLENDERTONE is a registered trade mark of BMR Ltd.
Bio-Medical Research Ltd., Parkmore Business Park West, Galway, Ireland.