Hæpi

HI-JULIA

Manuel d'utilisation
Instruction Manual
Manual de instrucciones
Bedienungsanleitung
Handleiding

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EN - 3-IN-1 ULTRASONIC TOOTHBRUSH
WITH WATER FLOSSER AND TARTAR
REMOVAL

Safety Instructions :

- **Toothbrushes can be used by children and by persons with reduced physical, sensory, or mental capabilities, or those lacking experience or knowledge, provided they are supervised or have been given instructions on how to**

use the device safely and understand the potential risks. Cleaning and user maintenance should not be carried out by children without supervision.

- **Water flossers can be used by children aged 8 and above and by persons with reduced physical, sensory, or mental capabilities, or those lacking experience or knowledge, provided they are supervised.**
- **Water flossers should only be used with cold or warm water, or with specific solutions defined in the manufacturer's instructions.**
- **Children must not play with the appliance.**
- **WARNING: Only use a Safety Extra Low Voltage power supply with the following specifications:**

 **5V == 1A / 5W**

*** Removable power supply part**

- **The device must be charged at a temperature of $20^{\circ}\text{C} \pm 5^{\circ}\text{C}$**
- **Do not expose the device or battery to excessive temperatures.**
- **Beware of the risk of short-circuiting the terminals of battery-powered devices or the battery itself with metallic objects.**
- **This device contains batteries that cannot be replaced. When the battery reaches the end of its life, the device must be disposed of appropriately.**

DISPOSAL AND RECYCLING:

WEEE Explanation: This symbol indicates that the product must not be disposed of with other household waste throughout the European Union. Improper disposal may harm the environment or human health, so please recycle it responsibly. This helps support the sustainable reuse of material resources. To return your used device, please use a return and collection system or contact the retailer where the product was purchased. They can dispose of the product in an environmentally responsible way.

Precautions for use

- Do not aim the water jet into the nose, ears, or eyes.
- Do not immerse the handle in water.
- Do not use if the device is damaged or wet.
- If you have dental or oral health conditions, consult a healthcare professional.

Additional Safety Warnings:

1. Never place the USB-C cable in water or any other liquid. Do not store it near bathtubs or sinks where it could fall. Do not touch it if it has fallen into water.
2. Avoid bending, pinching, or crushing the USB-C cable.
3. Do not use in the bath or shower.
4. This device contains a non-replaceable battery. There are no user-serviceable parts. Do not open the device.
5. This toothbrush is a personal hygiene device. It should not be used by multiple individuals in a clinic or institutional setting. For hygiene reasons, each user must have their own base and brush head.
6. Small parts may detach – keep out of reach of children.
7. To prevent brush head breakage, which may pose a choking hazard or damage teeth:
 - o Before each use, check that the brush head is securely attached. Do not use if it is loose.
 - o Replace the brush head every 3 months or sooner if worn.
 - o Properly clean the head after each use.

WARNING – To reduce the risk of burns, electric shock, fire, or injury:

- 1.**Proper use:** Follow the instructions. Use only accessories recommended by the manufacturer.
- 2.**Cord caution:** Keep the cord away from hot surfaces.
- 3.**Do not drop or insert objects:** Do not drop the device or insert anything into its openings.
- 4.**Check condition:** Do not use the device if the cord or plug is damaged, if it malfunctions, has been immersed in water, or is otherwise damaged. In this case, contact customer service.
- 5.**Not for outdoor** use or in environments with aerosol sprays or oxygen tanks.

1- Product Overview

Julia is a sonic toothbrush with integrated functions:

- Sonic brushing
- Interdental water jet
- Scaling head

It ensures complete oral hygiene in a single gesture.

2- Box Contents

- 1 Julia toothbrush
- 2 brushing heads
- 1 water flosser
- 1 tartar removal
- 1 water tube
- 1 travel case
- 1 wall mount
- 1 quick start guide
- 1 user manual
- 1 USB-C cable

3- Operating Modes

Fonction	Available modes	Niveaux de puissance
BRUSH	Clean/White	Normal (42 000) Strong (47 000)
JET	Continuous/ intermittent	Normal (60-85 PSI) / Strong (80-100 PSI)
COMBO	Brushing + Flossing	Possible with brush head attached and tube connected

- Easy combination via a single button
- Child-safe mode

4-Charging and Battery

- Fast charging via USB-C
- 3.7V 1100 mAh lithium battery
- Battery life: 15 days
- Charging time: 5-6 hours
- LED indicators for: mode, speed, battery

5- First Use

1. Fully charge Julia for at least 3 hours using the USB-C cable provided in the box. Ensure the device is properly connected to a stable power source. The LED indicator will show the charging level.
2. When all three bars are lit, the device is fully charged.
3. Attach the desired head:
 - Brush for standard cleaning
 - Jet for interdental cleaning
 - Scaling head for localized tartar
4. Fill the cup with water and connect the tube to activate the jet function.

Mode Setting:

To change the mode, press **brush and floss**

To change the speed, press **speed**

The toothbrush will automatically start in the last mode used.

Brush Heads:

When used properly – i.e. twice a day for two minutes – we recommend replacing the brush head every 3 months.

If the bristles flatten before 3 months, it means you're applying too much pressure to your teeth and gums while brushing.

6- Daily Use

- **Brushing:** Brush for 2 minutes with auto-timer (beep every 30 seconds)
- **Water flossing:** Aim the jet between the teeth, keep your lips slightly closed, and let the water flow into the sink
- **Tartar removal:** Use the special head to remove

visible tartar (once a week)

7. Cleaning and Maintenance

After each brushing, rinse the brush head under running water while keeping the handle on. Remove the brush head and clean both the handle and brush head separately. Ensure that water flows through the handle while it is on, as well as through the brush head. Wipe all parts dry and let them air dry before reassembling your toothbrush.

When the rechargeable battery is low, the indicator light will show only one bar. When the battery is fully discharged, the motor will stop, and the toothbrush must be charged for at least 5 minutes before it can be used again.

Note: If the battery is completely discharged (screen doesn't turn on), charge for at least 30 minutes.

Note: Store the toothbrush at room temperature for optimal battery life.

Warning: Do not expose the toothbrush to temperatures above 50°C (122°F).

Replacement Brush Heads Warranty

Haepi's warranty will be voided if damage to the rechargeable electric toothbrush handle is caused by the use of replacement brush heads from brands other than Haepi. Haepi does not recommend using third-party brush heads.

- Haepi has no control over the quality of non-Haepi brush heads. Therefore, we cannot guarantee their cleaning performance as originally promised when purchasing the Haepi rechargeable handle.
- Haepi cannot ensure proper fit of third-party brush heads.
- Haepi cannot predict long-term effects of using third-party brush heads on handle wear.

All Haepi replacement brush heads carry the Haepi logo and meet Haepi's high-quality standards. Haepi does not sell brush heads or spare parts under any other brand name.

8. Technical Specifications

- Materials: PC, ABS
- Dimensions: 36 x 36 x 185 mm
- Weight: 320 g
- Water pressure: 60 to 100 PSI

- Water jet frequency: 1400–1600 pulses/min
- Brushing frequency: 42,000 to 47,000 movements/min
- Power supply: 5V $\equiv$ 1A / 5W

9. Warranty & Support

- 2-year limited warranty
- Contact: info@distec.be

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