ID **Italian**
design

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EN PRECAUTIONS

- Please read the instructions carefully before using the appliance. To reduce any risk of accident, take the following precautions:
1. Always unplug the appliance immediately after use.
 2. Never use the appliance while you are in the bath.
 3. IEC 60417 - 5582: Do not place or store the appliance anywhere where it could fall or be pushed towards the bath or sink.

4. Do not immerse or drop into the water or any other liquid.
5. If the appliance falls into the water, unplug it immediately. Do not place your hand in the water.
6. Keep the appliance away from water.
7. Do not handle the appliance with wet hands.
8. An appliance should never be left unattended while it is plugged into the mains.
9. This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children must be watched to make sure that they do not play with the appliance.
10. Use this appliance only for its intended purpose, as indicated in this manual. Do not use any accessory that is not recommended by the manufacturer.
11. Do not use the appliance if you are drowsy.
12. Do not drop the appliance or insert any object in the openings or tube.
13. Disconnect the instrument with power supply before cleaning it.
14. Clean the residue of gels with a soft damp cloth or paper towel immediately after each use.
15. Do not use abrasive products for cleaning.
16. This care instruments suitable for private use only, do not share with your family or friends.
17. Preserve the device in a cool, dry place.
18. If the appliance, cord, or parts appear damaged or are not working properly, do not use or attempt repair; it should only be repaired by a qualified technician.
19. Fully charge the battery before first use.
20. Make sure that the voltage of the power supply corresponds to that of the appliance.
21. Do not pull on the cord and make sure it does not get tangled.
22. Keep cord away from hot and sharp surfaces.
23. If your skin is under nursing, or under any unhealthy condition, please consult your doctor before using the device.
24. IMPORTANT: do not use on the face.
25. For people with loose muscles, this device can leave a bruised mark after use.
26. In case of bruising occurs, stop to use it immediately. Normally, it will gradually disappear after 1 week.
27. Do not use on varicose veins.
28. Do not use during pregnancy or menstruation.
29. Do not use for more than 25 minutes per session.
30. Let the skin rest for a few days before the next session.

EN HOW TO USE

- Charge the device using the supplied USB cable.
- Clean the skin and apply plenty of body oil or lotion on it to help the device slide over the skin.
- Before each use, insert a round filter into the massage head with the help of tweezers and put the cap on. Place the second square filter in the filter holder of the body box to prevent lotion from entering the massager.
- Press the power button; select with the left arrow the heat level and with the right arrow the desired suction level.
- **IMPORTANT NOTE:** start the massage with the lowest suction level and increase it, little by little, according to your needs; the magnetotherapy function is continued.
- Slide the massager over your skin making slow and continuous movements from top to bottom on legs, buttocks, waist, abdomen and arms.
- Once the massage is finished, turn off the device: it has an automatic shutdown after 25 min of use.
- Remove the filters and clean them with warm water after each use; let them air dry before using them again.

EN TECHNICAL CHARACTERISTICS

- INPUT: DC 5V/2A.
- BATTERY: 3.7V / 1800mAh / 15W.

EN ENVIRONMENT

Disposal of old Electrical & Electronic Equipment (Applicable throughout the European Union and other European countries with separate collection programs. This symbol, found on your product or on its packaging, indicates that this product should not be treated as household waste when you wish to dispose of it. Instead, it should be handed over to an applicable collection point for the recycling of electrical and electronic equipment. By ensuring this product is disposed of correctly, you will help prevent negative consequences to the environment and human health, which could otherwise be caused by inappropriate disposal of this product. The recycling of materials will help to conserve natural resources. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

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Disposal of old Electrical & Electronic Equipment (Applicable throughout the European Union and other European countries with separate collection programs). This symbol, found on your product or on its packaging, indicates that this product should not be treated as

household waste when you wish to dispose of it. Instead, it should be handed over to an applicable collection point for the recycling of electrical and electronic equipment. By ensuring this product is disposed of correctly, you will help prevent negative consequences

to the environment and human health, which could otherwise be caused by inappropriate disposal of this product. The recycling of materials will help to conserve natural resources. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

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