ID Italian
design

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(EN) PRECAUCIONES

Please read the instructions carefully before proceeding with the use of this device. To reduce the risk of accidents, take the following precautions:

1. Do not use while taking a bath.
2. Warning: The IEC 60417 - 5582 symbol indicates that this appliance should not be used near a bathtub, shower, washbasin or other receptacle containing water.
3. Do not place or drop into water or other liquids.
4. If the appliance falls into water, unplug it immediately. Do not reach into the water.
5. Keep the device away from water.
6. The device must never be left unattended while it is plugged in.
7. Use this appliance only for its intended use as described in this manual. Do not use accessories not recommended by the manufacturer.
8. Never operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or otherwise damaged in any way, or dropped into water. Take the appliance to the nearest authorized service center for examination or repair.
9. Never block the openings of the appliance.
10. Keep the cord away from heat or hot surfaces. Do not wrap the cord around the appliance.
11. Never use if you are drowsy or half asleep.
12. Do not use outdoors or in places where aerosol (spray) products are being used or where oxygen is being administered.
13. This unit becomes hot during use. Do not allow hot surfaces to touch the skin.
14. Do not place the unit directly on any surface while it is still hot or plugged in.
15. This appliance can be used by children aged from 8 years and above and individuals with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
16. Do not tamper with the appliance. Opening the appliance will expose you to dangerous electrical components and pulsed light energy, which may cause serious physical injury and/or irreversible eye damage.
17. Clean the appliance only with a cloth slightly dampened with water, do not use corrosive products.
18. The IPL hair removal appliance is not as effective on white or light-coloured hair.
19. Do not use IPL hair removal on naturally dark skin.
20. Do not use the IPL hair removal device on tanned skin or after recent exposure to sun or UV equipment during the previous 28 days. Avoid exposure to sunlight for 4 weeks prior to IPL Hair Removal treatment and for 2 weeks after treatment.
21. Do not use IPL Hair Removal on white, grey or blond hair.
22. Do not use the IPL Hair Removal system on the scalp, eyebrows, chin, near the eyes, ears, nose, nipples, genitals, anus, as well as on areas with various marks, freckles, moles, warts, bruises, rashes, tattoos or permanent make-up, etc.
23. Do not use if pregnant or breastfeeding.

24. Do not use if you suffer from eczema, psoriasis, lesions, open wounds or evolving infection in the area to be treated. Wait until the affected area has healed before using the IPL hair removal device.
25. Do not use if you have a history of keloids, or if you have any pathology linked to photosensitivity, such as porphyria, polymorphous light eruption, solar urticaria, lupus, etc. or known sensitivity to light (photosensitivity) or caused by drugs, including non-steroidal anti-inflammatory drugs (e.g. aspirin, ibuprofen, etc.). Aspirin, ibuprofen, paracetamol, tricyclic, phenothiazine, thiazide, diuretics, sulfonamides, sulfonylureas, DTIC, fluorouracil, vinblastine, griseofulvin, alpha hydroxy acids (AHAs), beta hydroxy acids (BHAs), Retin-A R, Accutane R or topical retinoids. Or if you have recently been treated or are being treated with Alpha Hydroxy Acids (AHA), Beta Hydroxy Acids (BHA), Retin-A R, topical retinoids or azelaic acid.
26. Do not use if you have an abnormal skin condition caused, for example, by diabetes or other systemic or metabolic disease.
27. Do not use if you have been treated with Accutane R (isotretinoin) during the previous 6 months.
28. Do not use if you have been on a corticosteroid regimen for the previous 3 months.
29. Do not use if you have signs of pre-herpetic growth in the area to be treated, unless you have consulted a doctor and have been given preventive treatment prior to using IPL hair removal.
30. Do not use if you suffer from epilepsy.
31. Do not use if you have an active implant such as a pacemaker, urinary incontinence implant, insulin pump, etc.
32. Do not use if you have a history of skin cancer or areas subject to dermal neoplasia.
33. Do not use if you have undergone radiation or chemotherapy within the last 3 months.
34. Do not use if you have any other disease that your doctor considers to be a risk for the treatment.
35. Wait until the affected area has healed before using the IPL hair removal device.
36. If you are in any doubt about the risks to you of using the IPL hair removal device, consult your GP or dermatologist.
37. Avoid side effects! Do not treat the same area of skin more than once per session. Avoid using a pulse on the same area.
38. Never look directly into the light emitted by the flash lamp. Wear the included safety goggles whenever using the IPL epilator.
39. Never point the application surface of the IPL epilator to emit a light pulse in an open space.
40. Always ensure that you point the treatment surface towards the skin and that it is in full contact with the skin before delivering a pulse.
41. Caution, do not touch the lens of the applicator surface with your finger as it becomes hot.
42. Never use flammable liquids such as alcohol (including perfumes, deodorants or any other preparation containing alcohol) or acetone to clean your skin before using the epilator. It may cause a temporary change in pigmentation.
43. Avoid sun exposure for 4 weeks before an IPL hair removal treatment and again for 2 weeks after treatment.
44. After treatment, never expose yourself to unprotected direct sunlight for about fifteen minutes or continuous exposure to unprotected indirect sunlight for one hour.
45. Use adequate protection, such as clothing that provides sufficient coverage or a sunscreen (with a protection factor greater than 50). It should be used regularly 3 weeks before and after treatment.

46. It is important that you follow these instructions carefully if you are exposed to UV rays in a sunbed.
47. If a self-tanning cream has been applied to the skin, wait until the tanning effect of the cream has completely disappeared before using the appliance.
48. Always shave the area to be treated and make sure your skin is clean before using the IPL hair removal appliance.
49. Be sure to cover birthmarks and tattoos during IPL hair removal treatment and also cover any brown or black marks, freckles, nevi, warts or growths before using the IPL epilator.
50. If burns or blisters appear, STOP IMMEDIATELY.
51. If you are in any doubt about the risks to you of using the IPL hair removal device, consult your GP or dermatologist.
52. Carefully follow the instructions in this manual when selecting the power intensity during each session with the IPL epilator. Selecting an inappropriate power setting may cause side effects.
53. Do not use the IPL epilator on children.
54. **WARNING:** The plastic bags used to wrap this appliance or the packaging may be dangerous. Keep these bags out of the reach of babies and children. Do not use in cots, baby beds, prams or playpens. The thin film will stick to the nose and mouth and prevent breathing. A bag is not a toy.
55. To ensure additional protection, we recommend that you install a residual current device (RCD) with a rated tripping current of no more than 30 mA in the electrical circuit that supplies power to the bathroom. Ask your electrical installer for advice.
56. Stop the use of the appliance if you notice any defects or if it appears defective or emits an abnormal noise or odour. In particular, stop using the appliance immediately if the flashing lamp is cracked or stained.
57. Stop using the device when the ON/OFF/READY indicator light turns red and disconnect the IPL hair removal device adapter from the power outlet.
58. This appliance should not be used near a flammable mixture of anesthetic and air or with oxygen or nitrous oxide.
59. Always unplug the IPL hair removal device adapter from the power outlet immediately after use.
60. Unplug the IPL hair removal device before cleaning.

(EN) POSSIBLE SIDE EFFECTS

When your pulsed light epilator is used according to these instructions, side effects and complications are rare. However, all cosmetic treatments, even those designed for home use, involve some ups and downs, so it is important that you understand and accept these risks and possible complications when using a pulsed light hair removal system designed for home use.

If you overlap pulses or administer multiple pulses in the same patch:

A. Slightly uncomfortable skin sensation:

Although IPL hair removal is generally well tolerated, most people who use it experience a slight uncomfortable sensation during use, often described as a prickling sensation in the treated area. This sensation usually lasts throughout the application. Even after a few minutes. Any discomfort beyond this uncomfortable sensation is abnormal and means that you should either discontinue the treatment because you cannot tolerate it, or reduce the selected potency level.

B. Redness: Your skin may become red just after using the IPL epilator or within 24 hours after the treatment and the redness usually disappears within 24 hours. Consult your doctor if it persists for more than 2-3 days.

- C. Skin hypersensitivity:** The skin is more sensitive in the area that has been treated and may become dry and flaky.
- D. Burns or sores:** In very rare cases the skin may be injured or burned after application. It may take several weeks for complete healing and, in exceptional cases, a scar may remain.
- E. Scars:** Although these cases are very rare, permanent scarring may occur. Scars usually take the form of a superficial whitish lesion (hypotrophic). However, the scar may be large and red (hypertrophic) or large and extend beyond the injured area (keloid). More advanced cosmetic treatment may be necessary to improve the appearance.
- F. Pigmentation change:** The IPL hair removal device acts on the hair roots, in particular on the pigmented cells of the hair bulb. However, there is a risk of temporary hyperpigmentation (increased brown pigmentation or discolouration) or hypopigmentation (whitening) of the surrounding skin. The risk of pigmentation change is higher among darker-skinned people. Pigmentation change and discolouration are usually temporary. Permanent hyperpigmentation or hypopigmentation is rare.
- G. Excessive redness and swelling:** In rare cases, the treated skin may become very red and swollen. This reaction is more common in sensitive areas of the body. The redness and swelling should disappear within 2 to 7 days, during which time you should apply ice to the affected area regularly. Gentle cleansing is not a problem, but you should avoid any exposure to the sun.
- H. Infection:** A skin infection is extremely rare, but can occur after a burn or skin injury caused by IPL hair removal.
- I. Bruising:** In very rare cases, the use of IPL hair removal can cause some blue or purple bruising that may last for 5 to 10 days. As the blue fades, the skin may turn a rusty colour (hyperpigmentation) which could be permanent.
- Consult your doctor if you experience persistent discomfort or irritation during or after IPL hair removal treatment.**
- NOTE: If burns or blisters appear, STOP IMMEDIATELY.**

(EN) BEFORE USE

- Every hair on our body has three growth phases: anagen, catagen and telogen. These three phases play an important role in the process of the IPL hair removal device. The anagen phase is the hair growth phase, while the catagen and telogen phases are resting phases. This is the time it takes for the hair to nest and move from one part of the body to another. However, we can establish an average of between 18 and 24 months.
- Most hair follicles are in one of the two resting phases at any given time. Consequently, IPL hair removal will have no effect when they are in these resting phases. However, it will affect hair that is in the anagen phase. It is therefore important to understand the length of a complete follicular cycle in order to achieve long-lasting results with IPL hair removal.
- The amount of pigment in a person's skin determines not only the colour, but also the level of risk to which the user is exposed when using a form of hair removal that uses light, regardless of what it is. The darker the skin, the greater the risk.
- Select the intensity level with CAUTION:
 - The power intensity level designates the power of the pulsed light delivered to your skin during treatment, from the lowest(I) to the highest(5).

- The power indicator lights indicate the selected power level. Increasing the power level will improve the results, but will also increase the risk of side effects (
- **IMPORTANT:** Perform an intensity control test on a small area of skin 48 hours before use. If you do not notice any abnormal reaction or discomfort within 48 hours of the test, you can select the power level of your appliance for the first session using the table below:
- **SKIN PREPARATION:**
 - Your skin should be shaved, clean, dry and free of dust and deodorant residue.
 - Do not epilate (either with an electric epilator or wax) before using IPL hair removal. Epilation (removal of the hair by the root) will render the pulsed treatment ineffective.
- **SKIN TYPES: PHOTOTYPES (according to the fitzpatrick scale):**
 - I: Very white, milky skin- Never tans and always burns.
 - II: Light – coloured skin – Tans lightly and usually burns.
 - III: Beige skin – Tans moderately and sometimes burns.
 - IV: Medium brown skin – Tans easily and seldom burns.
 - V: Dark brown skin – Tans very easily and rarely burns.
 - VI: Black skin.
- For the first session, people with skin phototype I or II should set the device to level 3 maximum. Gradually, in subsequent sessions, if you do not experience any specific discomfort or discomfort, the power level can be increased to level 4 and then to level 5.
- In the first session, people with skin phototype III or IV should set the device to level 2 at most. Gradually, during the following sessions, if you do not experience any specific discomfort or discomfort, the power level can be increased to level 3 and then to level 4.
- In the first session, people with skin of skin type V should set the device to level 1 at most. Gradually, in subsequent sessions, if no specific discomfort or discomfort is experienced, the power level can be increased to level 2 and then to level 3.
- **The IPL hair removal device is not suitable for people with skin phototype VI.**

PHOTO TYPES	I	II	III	IV	V
SESSIONS Test (taken 48 hours before the first session)	Level 1	Level 1	Level 1	Level 1	Level 1
First session	Level 3	Level 3	Level 4	Level 5	Level 1
Subsequent sessions	Level 4,5	Level 4,5	Level 3,4	Level 3,4	Level 2,3

- Note:** The level of power intensity required may vary depending on the difference between a person's skin colour and the area to be treated. The skin on the arms, which is usually exposed to the sun more frequently, is for example darker than the skin under the armpits. Be careful to increase the power level for the lighter areas of your body. You will get better results, but you will also increase the risk of discomfort or even side effects.
- The first 4 sessions with IPL Hair Removal should be separated by 2 weeks.
 - Subsequent sessions should be spaced 4 weeks apart until the desired results are achieved.

- After that you will only need to use IPL hair removal when the hair grows back until the hair removal is long-lasting, leaving at least 4 weeks between sessions.

- **Note:** Do not treat the same patch more than once in two weeks, as this will not improve results and risks causing unwanted effects.

- **Note:** Due to hormonal factors and physiological changes, dormant follicles may become active. Therefore, maintenance sessions may sometimes be necessary.

- The first results are visible after 4 sessions and a very significant reduction in the amount of hair is achieved after 9 to 10 sessions. The result will vary from person to person, depending on the type of hair and the biological factors that determine the hair growth cycle. Some people achieve results more or less quickly than the average number of sessions we recommend. Lasting results cannot be expected after one or two sessions. Again, the length of the hair follicle resting cycle varies from one area of the body to another.

- The IPL epilator is designed with a built-in mechanism (contact sensor) that will protect your eyes. This was designed so that the pulse cannot be emitted if the epilator is pointed into an open space. The safety mechanism only allows the appliance to operate when the application surface is touching the skin. Do not look at the treated area when the appliance is flashing.

- IPL Hair Removal makes a slight fan noise, this is normal.

- Whenever you release a pulse, it is normal to feel some heat and pricking resulting from the light energy.

- The sensation of heat may persist for up to an hour after the epilation session.

- Right after use or within 24 hours after use, it is not uncommon for the skin to become slightly red or pink, especially around the hairs.

(EN) INSTRUCTIONS FOR USE

1. Take the appliance and the other IPL hair removal parts out of the box.
2. Plug the adapter cable into the device.
3. Plug the DC head of the adapter to the mains.
4. Press the button to switch on the appliance. The LCD display and the ready light will turn on and you will hear the sound of a fan. Your IPL epilator is now ready for use.

5. Select the power level that corresponds to your skin type (skin card). With the button, select the desired intensity level. To go up from level 1 to level 4, for example, press the button as many times as necessary, the appliance will go up from level 4 to level 5, then back from level 5 to level 4, from level 4 to level 3 and so on up to level 1.

6. WORKING MODES:

PRECISION/MANUAL mode: Press the release button each time you want to deliver a laser pulse. It allows safe and precise treatment of the less accessible and more sensitive areas of your body (face, underarms, bikini line, etc.).

GLISS/AUTOMATIC mode: Press and hold the release button to automatically release the laser pulses (continuous release of 36 pulses per minute). In this way, you can easily glide the device over the skin without interruption, in a slow, continuous movement, without lifting or lowering it. Please note that for safety reasons this continuous firing mode will stop as soon as you pull the sensor away from the skin.

7. Press the application surface firmly against your skin for full contact. When the application surface is in full contact with the skin. The READY indicator light will begin to flash.

Note: The appliance has a mechanism to protect your eyes, a contact sensor that prevents the appliance from pulsing if its application surface is not in full contact with your skin. Make sure, therefore, that the application surface is correctly applied to your skin.

8. Press the release mechanism. The device will release a pulse of light. You may feel a slight warmth and tingling sensation.

9. Remove the application surface from the treated skin area.

10. Move the application surface to another area to be treated. The pressure marks left on the skin by the application surface will help you to correctly align the epilator to deliver the next pulse, avoiding two pulses on the same area of skin or too much space between the areas (see section 'Covering the area to be treated'). Your IPL device will recharge immediately after each pulse. After 1.6 seconds (on level 1), the READY light will come on again. The device is now ready to deliver another pulse.

Warning: Never treat the same area of skin more than once per epilation session. Repeated treatment of the same area increases the risk of side effects.

Note: If you lift the appliance while holding down the release mechanism, the Gliss mode remains active. But for your safety, no impulse will be released until the appliance is in full contact with your skin. Make sure that the contact sensor that forms the contour of the applicator surface of the appliance is perfectly and firmly pressed against your skin to ensure continuous flashing operation in Gliss mode. If you break the contact between the applicator surface and your skin, you have 5 seconds to adjust the appliance to ensure full contact with your skin while pressing and holding the release button. If, after 5 seconds, the contact is not full, the Gliss mode will be deactivated and you will have to release the release button and press it again to restart the process.

Note: If you lift the device and release the release mechanism, the Gliss mode will be deactivated.

Note: Gliss Mode is suitable for people who already have some experience and skill with the device. With 36 beats per minute (on level 1), the Gliss Mode allows you to treat larger areas of the body more safely and quickly: legs, back, arms, torso, etc. For sensitive areas and the face, use the Precision Mode.

Note: If the appliance is not used for more than 10 minutes, it will automatically switch off. The "ON/OFF/READY" indicator light will illuminate orange. To use the appliance again, press the button and reset the desired power level.

11. Pulses should be delivered in rows, from one end to the other and proceeding systematically in successive rows. This method ensures better control of the pulses and will help you to better target the area in question and avoid treating the same patch twice or delivering pulses in the same area or delivering pulses too close together in the same area. The pressure marks left on the skin by the application surface will also help you to align the Epilator correctly.

IMPORTANT: STOP treatment IMMEDIATELY if your skin burns or blisters.

Area	Average duration	Mode
Half a leg	8 to 10 minutes	Gliss
Leg	15 minutes	Gliss
Bikini line	2 minutes	Precision
Underarms	1 minute	Precision

FACIAL DEPILETION:

1. It is essential not to look directly at the treated area during the flash and to wear protective goggles when IPL hair removal is used to treat the face.
2. Do not use for plucking eyelashes, head hair and eyebrows.
3. This product is not suitable for full beard depilation of the male beard.
4. For your comfort and safety, we recommend that you have someone assist you during the first few sessions when using the IPL epilator on your face.
5. The PRECISION mode is required to remove facial hair.
6. **DO NOT USE GLISS MODE TO TREAT THE FACE OR SENSITIVE AREAS.**
7. To treat the face, it is essential to wait 5 seconds between flashes.
8. The duration of the IPL hair removal session on the face should not exceed 10 minutes.

SKIN CARE AFTER THE DEPILETION SESSION:

1. Avoid exposing treated areas to the sun. Protect your skin with sunscreen throughout the treatment and for at least two weeks after the last session with your IPL hair removal device.
2. Do not use perfumed creams or deodorants on the treated areas after the sessions.

(EN) CLEANING AND MAINTENANCE

1. At the end of the IPL hair removal session, switch off the appliance with the off button and unplug the adapter from the mains.
2. After each session, gently wipe the surface of the device and, in particular, the application surface with a clean, dry cloth and a cleaning product formulated for cleaning electronic equipment.
3. Store in a dry and dust-free place.

(EN) TROUBLESHOOTING

1. My IPL epilator does not turn on.

Make sure that the adapter is correctly connected to your IPL epilator and plugged into the mains.

2. There is no impulse when I press the release mechanism.

Make sure that the application surface is firmly in contact with your skin and that the "ON/OFF/READY" indicator light is flashing; this light ensures that the device is in full contact with your skin. For your safety, the release mechanism will only be able to release a pulse if the application surface is firmly pressed against your skin.

- Make sure the "ON/OFF/READY" light is on. If after 10 seconds the 'Ready' light is still off, switch off your IPL epilator and switch it on again by pressing the button twice. If the problem persists, contact your local customer service centre.

- The ON/OFF/READY light is read and flashing. The appliance has gone into safety mode because it has detected that it is overheated or there is a critical problem (broken fan, condenser problem, etc.). In these circumstances, switch off the appliance and let it cool down for 15 minutes before switching it on again. If the problem persists, contact the customer service in your area.

- The "ON/OFF/READY" indicator light flashes green and red alternately. The adapter used to connect the IPL hair removal appliance to the mains is not the one supplied with the appliance or is defective. If the problem persists, contact your local customer service centre.

- The "ON/OFF/READY" indicator light is red. The appliance has reached the end of its service life. Do not attempt to open or repair the

appliance yourself, as this may not only invalidate the warranty, but may also cause serious physical injury and/or irreversible eye damage. Only certified repair centres are authorised to carry out repairs.

(EN) SPECIFICATIONS

Point of impact: 3.0 x 1.0cm=3.0cm².

Repetition rate:

- At level 1, one pulse every 1.6 seconds.
- At level 5, one pulse every 3.5 seconds

Luminous intensity:

- At level 1: 2.0 J/cm
- At level 2: 2.5 J/cm
- At level 3: 2.8 J/cm
- At level 4: 3.2 J/cm
- At level 5: 3.5 J/cm

Wavelength: 550-1200nm.

Required electrical installation: AC 100-240V, 4A.

(EN) ENVIRONMENT

Disposal of old Electrical & Electronic Equipment (Applicable throughout the European Union and other European countries with separate collection programs.)

This symbol, found on your product or on its packaging, indicates that this product should not be treated as household waste when you wish to dispose of it. Instead, it should be handed over to an applicable collection point for the recycling of electrical and electronic equipment. By ensuring this product is disposed of correctly, you will help prevent negative consequences to the environment and human health, which could otherwise be caused by inappropriate disposal of this product. The recycling of materials will help to conserve natural resources. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

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