### MODES

**MODE 1**  
2 min

**MODE 2**  
2 min

**MODE 3**

**MODE**  
3 sec

**ON** **OFF**

1 x

2 x

3 x

100%

24h

**Deutsch**

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If the toothbrush handle is dropped, the brush head should be replaced before the next use even if no damage is visible.

Replace brush (head) every 3 months or sooner if brush head becomes worn.

When unplugging, always pull the power plug instead of the cord. Do not touch the power plug with wet hands. This can cause an electric shock.

For tooth brushing only.

Use this product only for its intended use as described in this manual. Do not use attachments which are not recommended by the manufacturer.

If stickers are provided, follow instructions on the sticker sheet. Ensure valve on back of the handle is not covered.

Periodically check the entire product/cord/accessories for damage. A damaged or non-functioning unit should no longer be used.

Do not modify or repair the product. This includes fire, electric shock or injury.

Never insert any object into any opening of the appliance.

This toothbrush has a waterproof handle, is electrically safe and designed for use in the bathroom.

Periodically check the operation of the Smart Pressure Sensor by pressing moderately on the brush head during use.

If you use a medical device or have an implanted medical device, consult with the medical device manufacturer or a physician, if you are using Oral-B IO toothbrush (accessory) contain NfFeB magnets and generate magnetic fields. If you use a medical device or have an implanted medical device, consult with the medical device manufacturer or a physician, if you are using Oral-B IO toothbrush.

**30 Day Risk Free Trial**  
Satisfaction guaranteed or your money back!  
Full details at [www.oralb.com](http://www.oralb.com) or call: (UK) 0800 731 1792 (IE) 1800 509 448

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**English**

Welcome to Oral-B! Before operating this toothbrush, please read these instructions and save this manual for future reference.

**Note:** Content may vary based on the model purchased.

**IMPORTANT**

Toothbrushes can be used by children and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge, if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.

If the product/cord is damaged, take it to an Oral-B Service Centre.

Only use charger provided with your appliance.

**Brush heads**

Oral-B offers you a variety of different Oral-B brush heads that only fit your Oral-B IO toothbrush handle and can be used for tooth-by-tooth precision cleaning.

**Cleaning guidelines**

After each brushing, rinse your brush head under running water with the toothbrush switched on. Remove the brush head and clean the toothbrush handle and brush head separately. Make sure water runs inside the top of the handle while the handle is switched on, as well as inside the brush head. Wipe all parts dry and let them air-dry before reassembling your toothbrush.

**Environmental Notice**

Product contains batteries and/or rechargeable electric waste. For environment protection do not dispose in household waste, but for recycling in designated waste collection points provided in your country.

**Guarantee**

We grant 2 years guarantee on the date of purchase.

Within the guarantee period we will eliminate, free of charge, any defects in the appliance resulting from faults in materials or workmanship, either by repairing or replacing the complete appliance as we may choose. Depending on availability replacement of a unit might result in a different colour or an equivalent model.

This guarantee extends to every country where this appliance is supplied by Braun / Oral-B or its appointed distributor. This guarantee does not cover damage due to improper use, normal wear or use, especially regarding the brush heads, as well as defects that have a negligible effect on the value or operation of the appliance.

The guarantee becomes void if repairs are undertaken by unauthorised persons and/or if original Braun / Oral-B parts are not used.

To obtain service within the guarantee period, hand in or send the complete appliance with your sales receipt to an authorised Braun Oral-B Service Centre.

**Warranty replacement brush heads**

Oral-B warranty will be voided if electric rechargeable toothbrush handle damage is found to be attributed to the use of non-Oral-B replacement brush heads.

Oral-B has no control over the quality of non-Oral-B replacement brush heads, therefore, we cannot ensure the cleaning performance of non-Oral-B replacement brush heads, as communicated with the electric rechargeable handle at time of initial purchase.

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