

Silk·expert Pro 5 PL5431

2 years of smooth skin. Visible results in 2 sessions.¹

BRAUN

Feature highlights



Smart SkinProtect sensor
auto-adjusts the power to your skin tone.



10 intensity levels
for the right power to you.



3 modes
Sensitive, Balanced and Power



Full body² treatment in 10 minutes



400,000 flashes
for 47 years of full body treatments.³



Flashes up to 2.08 times per second
for an effortless treatment and better skin coverage.



Glides smoothly for an easy session



Gentle on intimate
Dermatologically and gynecologically tested.



Built-in UV filter
No safety glasses needed.



Corded use for non-stop power
Requires no batteries or charging.



Designed in Germany

Areas of use




Legs



Arms



Underarms



Chest



Back



Face



Bikini area



Intimate area



What's in the box

SMART FLEX HEAD GOLD

STANDARD HEAD

PRECISION HEAD

MINI PRECISION HEAD

VENUS EXTRA SMOOTH RAZOR

POUCH

8700216937948

¹ Following recommended regimen. Individual results may vary. ² Legs (upper & lower), underarm, bikini, upper lip. ³ When following recommended regimen.

Features in detail



Smart SkinProtect sensor

Skin tone varies across the body. Our best sensor reads your skin tone 80x per second, automatically adjusting each flash to your different body areas for a safe treatment.



3 treatment modes

Choose the Sensitive Mode for a gentle treatment of sensitive areas, the Balanced Mode for a balance between quick results and a gentle session, or the Power Mode for quick results.



Our fastest IPL

Treat your full body in just 10 minutes.¹ Requires only 2 sessions per month for a total of 3 months. Then touch-up, if needed.



400,000 flashes

Braun IPL Silk·expert Pro 5 comes with 400,000 flashes that deliver an equivalent of 47 years of full body treatments.²



Gentle on intimate

Braun IPL devices are both dermatologically approved and gynecologically tested, making them suitable for use on intimate areas.³

Accessories in detail



Smart Flex head

Thanks to flexible wings, the Smart Flex head adapts to the curves of your body for optimal skin contact to avoid missed hairs. Ideal for legs, arms and tricky areas.



Standard head

Our head attachments are designed for different body parts. The Standard head is well suited for the treatment of underarms.



Precision head

Use the Precision head for a gentle treatment of intimate areas.



Mini Precision head

Designed to treat small areas like face and toes with precision.



Gillette Venus Razor

Before IPL treatment, remove all visible hair in the area you want to treat with the razor. This will help the light pulse to penetrate the hair root more easily.

¹ Legs (upper & lower), underarm, bikini, upper lip ² When following recommended regimen. ³ Not to be used by women on: scalp, eyes, forehead, nipples, on mucous membranes such as the labia minora, vagina, anus; and men on: scalp, face, neck; nipples, penile shaft, scrotum, anus.