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
Dear Customer,

Thank you for purchasing this product. To be completely satisfied with your product, we recommend reading this user manual carefully.

If you have any difficulties or need more information, you can find our help centre by scanning the QR code opposite with your smartphone. If you cannot find the answer to your question, you can also contact our After Sales Service: support.uk@terraillon.com

WHY CHOOSE SLIM CARE ?

The SLIM CARE belt combines EMS (Electrical Muscle Stimulation) and TENS (Transcutaneous Electrical Nerve Stimulation) with soothing heat for comfort and effective toning. Made from skin-friendly Lycra, the belt is ultra-thin, breathable, and fits discreetly under clothing. Its elastic design adapts to most waist sizes and can be worn on the abdominal or lumbar area.

Ideal for fitness support or workout enhancement, the belt uses medium and low-frequency waves to stimulate muscles and relieve tension. The integrated heating function adds extra relaxation, used alone or with electrostimulation programs.

This manual provides guidance on program modes, intensity, and duration, which can be adjusted for your comfort.

Electrical Muscle Stimulation (EMS)

Electrical Muscle Stimulation uses controlled electrical impulses to make muscles contract and then relax, just like during physical activity. These repeated muscle contractions help strengthen and tone muscles, improve body shape, and support slimming goals. By mimicking exercise, EMS offers a form of passive workout that can complement training or assist recovery. It also helps boost circulation, supporting muscle health and faster recovery.

Transcutaneous Electrical Nerve Stimulation (TENS)

TENS delivers mild electrical impulses through the skin to gently stimulate nerve endings. This process helps reduce muscle tension, soothe soreness, and promote deep relaxation by encouraging natural endorphin release. TENS is ideal for reducing stress in muscles, supporting recovery after physical activity, and enhancing overall comfort and well-being.

Soothing Heat

To enhance the electrostimulation effect, or as a standalone feature, the heat function provides soothing warmth that helps relax muscles and ease discomfort. Heat therapy also supports the body's natural recovery process after physical exercise, making it an ideal complement to your fitness routine.

Designed for Everyday Use

SLIM CARE features a unique slim profile and soft, comfortable fabric, making it easy to wear anytime and during most of your daily activities. Its modern yet discreet design allows you to wear it under your clothes or use it during physical exercise without compromising comfort or style.

CONTRAINDICATIONS

The following contraindications must be observed before using SLIM CARE, also called the device in this instruction manual. Read the operating instructions carefully and keep them for future reference. If you pass this device on to others, please also provide this manual.

SLIM CARE is intended for healthy adults between 18 and 60 years of age. If you have any concerns or health-related questions, consult your physician before use.

DO NOT USE THIS DEVICE IF YOU:

- Have a cardiac pacemaker, implanted defibrillator, or any other implanted metallic or electronic device.
- Are in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms).

- Are near instruments that emit shortwave or microwave.
- Have a history of prolonged drug use or injections.
- Have gastrointestinal bleeding, liver cirrhosis, or diseases affecting the spleen.
- Are experiencing hemorrhage during the acute or unstable phase. For any residual effects or sequelae, use only under medical supervision.
- Have acute cardiovascular or cerebrovascular conditions.
- Have open wounds, rashes, swollen, infected or inflamed areas, or cancerous lesions.
- Have skin allergies, heart disease, hypertension, malignant tumors, severe cerebrovascular diseases, cancer, epilepsy, or other functional disorders.
- Have a fever, as the heating function can increase body temperature.
- Are insensitive to heat, you must be careful when using the appliance, as it has a heating function.
- Plan to use the device on the head, chest, heart, or any body part not specified in this manual.

CONSULT YOUR DOCTOR BEFORE USE IF YOU:

- Have been pregnant or recently gave birth.
- Experience abdominal pain or have a serious illness or chronic disease.
- Have recently undergone surgery, suffer from diabetes, or have a muscle injury.
- Plan to use the device for health recovery purposes.

DO NOT USE ON THESE INDIVIDUALS:

- Pregnant women.
- Children or infants.

DO NOT USE DURING THESE ACTIVITIES:

- Bathing, showering or near water sources.
- Sleeping.
- Driving, operating heavy machinery or any activity where electrical stimulation could increase the risk of injury.

PRECAUTIONS

GENERAL INFORMATION

- Long-term effects of electrical stimulation are unknown.
- Effectiveness depends on the user.
- This should not be used together with other products.
- Keep out of reach of children (contains small parts that may cause choking or be swallowed).
- Prevent misuse of connecting wires to avoid strangulation accidents.

PRODUCT SAFETY:

- Do not modify, repair, or disassemble the device or remote without manufacturer authorization.
- Check the device and remote for damage before use.
- Clean the device with a dry, soft cloth (dust and lint may affect performance).

USAGE PRECAUTIONS:

- Stop using immediately if the device malfunctions or if you feel discomfort, dizziness, or palpitations.
- For first-time users, start with low stimulation intensity and adjust the mode to avoid discomfort.
- Do not move or change the device's position during use (risk of strong irritation).
- Do not use for more than 30 minutes per day for electrostimulation and heating functions.
- Remove any belly piercing or jewelry near usage zone before use.
- Use only the remote approved by the manufacturer. The use of unapproved remote may damage the device or injure users.
- If the device belt is wet, let it to dry before use.

POSSIBLE ADVERSE REACTIONS:

- This product has a heating function, and prolonged continuous use (more than 30 minutes) may cause low-temperature burns or skin irritation.


- Muscle fatigue may occur if recommended usage time is exceeded.
- Stop using and consult your physician if adverse reactions occur.

This device may only be used for the intended purpose described in this manual. The manufacturer cannot be held liable for incidental, consequential or special damage caused by misuse or abuse.

PACKAGE CONTENTS

Belt
x1

Remote control
x1

Extension
x1

USB-C cable
x1


Instruction manual
x1


Quick start guide
x1

PRODUCT DIAGRAM

Front


Back

Terraillon®

Remote Control

PREPARATION BEFORE USE

1. CHARGING THE PRODUCT

Before the first use, charge the device to maximum capacity.

Connect the USB-C cable to its port and connect the other end to an adapter.

When charging, the remote control is blinking in orange. When battery is fully charged (about 2 to 3 hours), the orange light becomes steady. USB cable can be disconnected and device operated. When the battery is low, the device emits three beeps.

Fully charged device can last about 3 to 4 days with 30 minutes usage per day. If the battery remains unused for a long time, its service life may be shortened. To maintain performance, it is recommended to fully charge the battery at least once a month.

2. SLIM CARE SET UP

Before you begin, make sure your abdomen or lower back is clean and slightly moist, this helps the electrodes work properly. If your skin is too dry, the connection might feel weaker.

Wrap the belt comfortably around your waist so the electrodes rest against your skin. Adjust the velcro strap for a snug fit, and use the extension if you need extra length. The electrodes can be placed on either your abdomen or lower back.

You can either select your programme before attaching the remote control to its location, or attach the remote control first and then choose your programme. If this is your first time, start with a low intensity and increase gradually as you get used to it.

You can wear the belt while going about your day or during light physical activity. It's designed to fit right into your routine.

3. OPERATING INSTRUCTION

1. **Turn on** the remote control by pressing the Power button . A short sound will confirm activation, and the indicator light will stay white. The device starts with default settings: intensity at 0, heating level at 1, and a timer set for 15 minutes.
2. **Select a mode** by pressing the Mode button  to cycle through options 1 to 6. A buzzer will indicate the selected mode:
 - 1 beep for Mode 1
 - 2 beeps for Mode 2
 - 3 beeps for Mode 3, and so on
3. **Adjust intensity** using the "+" or "-" buttons. There are 16 intensity levels available. Start low and increase gradually for comfort.
4. **Activate heating** by pressing the heating button  to cycle through: Low, Medium, High, and Off. Each press changes the level, and the buzzer will confirm:
 - 1 beep for Low : around 38°C
 - 2 beeps for Medium: Around 39°C to 41°C
 - 3 beeps for High: Around 42°C to 45°C
 - Long beep when turning heating off

The belt will automatically shut down after 15 minutes. To stop earlier, press and hold the Power button .

For best results, use once or twice daily, with a maximum of 30 minutes per day for EMS/TENS and heating combined.

To support your journey, we've curated a selection of standard programs. These are general guidelines, and you're free to customize or create your own based on your needs and comfort level. For first-time users, always start with low intensity.

You can adjust the intensity based on your experience with electrostimulation and your comfort level. The intensity levels are defined as follows for all programs, feel free to customized it based on your personal experience:

- Very low: 1 to 4
- Low: 5 to 7
- Medium: 8 to 11
- High: 12 to 16

Programs	Combinations of modes	Intensity	Duration/day
Beginner Program	Week 1/2: Mode 3	Very low → Low	15 min
	Week 3/4: Mode 5	Low	15–30 min
	Week 6: Mode 1	Low	30 min
Toning Program	Week 1/2: Mode 4	Low	15 min
	Week 3/4: Mode 3	Medium	30 min
	Week 5/6: Mode 2	Medium → High	30 min
Gentle Sport Program	Week 1/2: Mode 4	Low	30 min
	Week 3/4: Mode 5	Medium	
	Week 5/6: Mode 2	Medium	
Active Sport Program	Week 1/2: Mode 4	Medium	30 min
	Week 3/4: Mode 6	Medium	
	Week 5/6: Mode 2	High	
Warm-up Program	Week 1/2: Mode 4	Medium	15 min
	Week 3/4: Mode 5	High	
	Week 5/6: Mode 6	High	
Gentle Massage & Relaxation	Use modes 1, 3, or 5 separately, or alternate between them.	Low, with or without heating function.	30 min
Energetic Massage / Relaxation	Use modes 3, 6, or 2 separately, or alternate between them.	Low to Medium, with or without heating function.	30 min

If you use both fitness and relaxation programs on the same day, do not exceed 30 minutes, as this may cause adverse effects.

The heating function can be included with relaxation programs, and you can activate it based on your needs and comfort. EMS/TENS and heating can be used independently, allowing you to use the belt for up to 1 hour (30 minutes EMS/TENS and 30 minutes heating).

FREQUENTLY ASKED QUESTIONS

Can lying-in woman use this device?

Yes, she can use the device after delivery (around 6 weeks), but must consult a doctor first to assess her specific condition. In the case of a cesarean section, use is only allowed after obtaining the doctor's approval later (around 3 months) in the recovery period.

I'm experiencing uncomfortable irritation. How can I relieve it?

1. Check if the belt has been worn correctly and adjust electrodes placement if needed, with the belt off.
2. If necessary, apply small amount of water on the surface of patches to make you feel better. Please note, do switch off the power or do not splash water on the controller when operating.
3. If your skin is too dry, apply hot compress function for 5 to 10 minutes and then use electric pulse massage function.

My skin become red after using the device, is there any problem?

This is a normal reaction. Low-frequency stimulation increases blood flow, which can make the skin appear red. Additionally, tightening the waist belt may leave temporary marks on the skin. Don't worry, these marks will gradually fade.

If the marks are too pronounced, loosen the belt and select a lower intensity level. If the irritation persists, stop using it and consult a doctor.

Will using the SLIM CARE cause muscular aches?

Just like all sports, using the device may cause different degrees of aches on body muscles. It is a normal phenomenon and will last 1 to 2 day(s). Please adjust use frequency according to your own situations.

TROUBLESHOOTING

PROBLEM	POSSIBLE CAUSES	SOLUTION
The remote control won't power on	Internal rechargeable battery is low.	Charge the controller using a USB-C cable.
	Charging issue.	Ensure the charging source outputs is 5V; indicator light should blink when charging.
Unable to charge or no charging display	USB-C cable not properly connected.	Check both ends of the cable are securely inserted.
	Insufficient charging power.	Verify power source output is normal.
Normal display, but no response or sensation	Electrodes are not in direct contact with skin.	Ensure electrodes touch skin directly.
	Clothing is between electrodes and skin.	Remove clothing or place belt under clothes.
	Strength not adjusted.	The default setting is 0 strength. Adjust the strength slowly as sensation may depend by users. If there is still no feeling, please use two finger to touch two electrodes, if no sensation please contact maintenance.
Sense of piercing pain	Electrodes not in proper contact.	Wipe skin with damp towel and ensure electrodes fit snugly, thus adjust settings.
	Usage duration is too long.	Limit usage to maximum 30 minutes per day.
	Strength is too high.	Reduce strength to a comfortable level.
	Skin sensitivity.	For allergies, please consult your doctor for advice.
	Cold weather/dry skin.	Wipe skin and electrodes with damp towel.
	Used after staying up late.	Use the next day if electrolyte levels are low.
Skin irritation	Electrodes not in direct contact with the skin.	Ensure proper contact and adjust strength.
	Dirt or oil on electrodes.	Clean electrode surfaces.
	Low battery capacity.	Use external power source if needed.
Red dots appeared after use	Strong stinging sensation during use.	Stop using the device; symptoms usually disappear within 3–5 days.

MAINTENANCE & CARE

CLEANING

1. When the massage belt fabric becomes dirty, use a dry soft cloth or dry towel to clean it.
2. After each use, please clean the electrodes. You can use a soft cloth or towel slightly dampened with water to clean it. When the conductive piece is particularly dirty, you can use a soft cloth slightly dampened with medical alcohol (75% concentration) to clean or disinfect it.


STORAGE

1. Please do not store it in direct sunlight, high temperatures, or damp areas.
2. Please store it in a dry and well-ventilated place.
3. Please do not disassemble, repair, or modify this product without authorization, which may cause accidents or malfunctions, except for battery disposal as below.
4. Please store it in a place that children cannot reach.

DISPOSAL

To prevent environmental pollution, do not discard the SLIM CARE belt or its accessories with household waste. When they can no longer be used, remove the battery and dispose of them according to local environmental regulations.

- The built-in lithium batteries contains substances that may contaminate the environment. Before disposal of it, remove the built-in lithium batteries.
- Always switch off or disconnect the power before disassembling the remote controller and removing the lithium batteries.
- If the device is damaged or shows any leakage, do not attempt to open it. Please contact the company after sales service for assistance.
- To remove the lithium battery, proceed as shown in the diagram below, and disassemble the product as per the steps below.

TECHNICAL SPECIFICATIONS

Model	SLIM CARE
Pulse	Pulse repetition frequency: 1Hz -1000Hz; Pulse modulation frequency: 4KHz-10KHz Pulse Width: 20 μ s -400 μ s; Rated Load Impedance: 500 Ω , permissible deviation of $\pm 10\%$
Time	Up to 30 min
Heating	3 level: Within 38°C to 45°C (under the environment of 23°C $\pm$ 2°C)
Dimensions	Box: 29,8 x 13,3 x 4,9 cm Belt: 96,5 x 6/11,8 x 0,01 cm Belt extension: 33 cm x 6 cm remote control: 10 x 1,5 cm
Weight	About 150g
Power sources	Charging Source: 5V DC 1A (USB) Charging time: around 2-3 hours Input Power: 5W Lithium Battery of Host: 3.7VDC 2000mAh Battery life: around 500 times charge and discharge, about 3 years

Operating environment	Ambient Temperature: +5°C-+40°C; Relative Humidity: ≤80%; Atmospheric Pressure: 86kPa -106kPa
Storage & Transportation conditions	Ambient Temperature Range: -20°C-+55°C; Relative Humidity Range: < ≤93%; Atmospheric Pressure Range: 50 kPa -106 kPa

EXPLANATION OF SYMBOLS

The following symbols may appear in the instructions for use, on the packaging, or on the product:

Do not wash

Do not dry clean

Do not machine wash

Do not bleach

Do not iron

Read the instructions

This product meets the requirements of current European and national directives

Disposal in accordance with the European directive WEEE (Waste Electrical and Electronic Equipment)

WARRANTY AND ENVIRONMENTAL PROTECTION

This appliance is guaranteed against material and manufacturing defects for two years. During this period, the product will be replaced free of charge (proof of purchase must be presented if the product is under warranty). This warranty does not cover damage resulting from accidents, improper use or negligence.

If you have a complaint, please contact the store where you purchased your product first.

Waste electrical products should not be disposed of with household waste.

Please recycle where facilities exist. Check with your local authority or retailer for recycling advice.

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Terraillon SAS France & Headquarters
Terraillon Asia Pacific Ltd
Terraillon Corp USA