
SKIN-BOOST

EZICLEAN

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INSTRUCTIONS

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Product Overview

Our LED Beauty Mask is made with soft, skin-friendly liquid silicone on both the inner and outer layers, as well as the eye mask. It is designed to be odorless and features adjustable straps to ensure a comfortable fit for all face shapes. Equipped with 54 advanced 4-in-1 LED bulbs, the mask emits four wavelengths of light: near-infrared, red, blue, and yellow, delivering a total of 216 light beams. The energy levels and treatment duration can be customized using the included controller.

The innovative hollow design reduces the mask's weight, eliminating the discomfort of sagging or pressure, making it more comfortable for prolonged wear.

The suspended design, supported by eye pads, ensures the mask remains approximately 15mm away from your skin, allowing the LED light to distribute more evenly across your face. This feature resolves the constraints and stuffiness typically associated with traditional masks, enabling you to breathe comfortably during treatment.

With a one-piece molding process, the LED light source is seamlessly embedded into the liquid silicone, simplifying the design while effectively protecting internal components and extending the mask's lifespan.

The mask is IP67-rated for dust and water resistance, allowing it to be rinsed directly under running water for easy and hygienic cleaning.

Product Specifications

LED Therapy Mask

Controller

Type-C Charging Cable

Adjustable Strap

Model	MMZ01
Net Weight	93G
Dimensions	210*210*160 MM
Input	DC 5V $\overline{\text{---}}$ 1 A
Power	4W
Battery capacity	Li-ion 2600mAh

Benefits of LED Technology for the Skin

The mask uses different light wavelengths designed to target the skin's surface layers, helping to support its natural renewal process, enhance radiance, and promote the appearance of firmer, more resilient skin.


- **Infrared Light (850 nm) + Red Light (630 nm)** : Helps improve the appearance of skin firmness and elasticity by supporting collagen and elastin. It visibly smooths the look of fine lines and wrinkles, refines skin texture, and helps reduce the appearance of pigmentation spots.
- **Yellow Light (590 nm)** : Provides comfort to delicate skin, helps promote a more even-looking complexion, enhances natural radiance, improves the absorption of cosmetic skincare products, and helps reduce the appearance of surface redness.
- **Blue Light (415 nm)** : Targets blemish-prone skin, helps purify the complexion, minimizes the appearance of enlarged pores, and contributes to rebalancing combination to oily skin.

Your Daily Skincare Ally

- Reduce the appearance of blemishes and breakout-prone areas
- Visibly smooth fine lines and wrinkles
- Soothe and comfort the skin
- Improve the appearance of skin firmness and elasticity
- Even, brighten, and enhance the complexion
- Support the skin's natural regeneration process

How to Use

1. Controller Overview

2. Charging the Controller


- Before your first use, **fully** charge the controller.
- During charging, the indicator light will blink. Once all lights are solid and no longer flashing, the controller is fully charged.
- After charging, disconnect the controller.
- Do not use the mask while charging.
- Only use the provided USB cable to avoid malfunction.

The mask does not have a built-in battery: the battery is integrated directly into the remote control.

To use the device, first charge the remote control via its USB-C port. Then, connect the remote control to the mask using the same cable, plugging the USB end into the remote control and the USB-C end into the mask.

3、Wearing the Mask


①

Velcro tie

Ensure the outer surface of the mask is facing you. Thread the strap with Velcro through the left ear hole and secure it.

②


Thread the other side of the strap through the mask, adjust to fit comfortably, and secure it.

③

Wear the mask from the top of your head, adjusting the nose bridge for comfort. If you have a ponytail, thread it through the strap hole.

④


Adjust the strap as needed to ensure the mask fits snugly but comfortably.

Always wear the mask before switching on the LEDs to prevent direct exposure of the eyes to the light.


4. Connecting the Mask to the Controller

Attach the mask to the controller using the provided cable.


5、Powering On/Off

Press and hold the power button for 1.5 seconds to turn on the device. The display screen will light up.

6、Selecting a Mode

Press the “M” button to switch between light therapy modes: red, blue, or yellow light. Once a mode is selected, the timer will start automatically.

Near-infrared + Red Light
8 50nm+ 630nm

Blue Light
41 5nm


Yellow Light
590nm

7. Adjusting Energy Levels

Short press the power button to cycle through three energy levels. The higher the level, the brighter the light and the stronger the treatment.

Recommendation:

Begin with the lowest level for your first use. If your skin is not sensitive, gradually increase to higher levels. Sensitive skin should use the lowest setting.


8. Setting Treatment Duration

Press and hold the “M” button for 1.5 seconds to set the treatment duration. Choose from 10, 15, 20, or 30 minutes.

We recommend using the mask for 10 minutes daily, not exceeding a total of 30 minutes per day.

You may also combine treatments, such as 10 minutes of red light followed by 10 minutes of blue or yellow light.

Steps for Safe Use

1. Sensitivity Test

Before your first treatment, perform a sensitivity test:

- Place the device on your inner arm.
- Turn on the device and run a 10-minute treatment.
- If you experience discomfort, stop immediately.
- After treatment, wait at least 6 hours. If redness or itching occurs, do not use the mask on your face.

2. Starting Your Treatment

- Cleanse your face with a gentle cleanser and use the mask on clean, dry skin without applying any creams or serums, to ensure optimal LED light penetration.
- Connect the mask to the controller.
- Press and hold the power button for 2 seconds to turn on the device, select your mode, and begin your treatment.
- After treatment, follow your usual skincare routine.
- Recommended use: 3-5 times per week, starting with every other day.

Cleaning and Maintenance

- Clean the mask before and after each use with a soft, dry cloth or wet wipe.
- The mask can be rinsed under running water, but ensure the controller is unplugged and disconnected.
- Do not rub or twist the mask forcefully.
- Store the device in a dry environment, away from direct sunlight.

RECOMMENDED USE

- Start with 10-minute sessions.
- Use 3 to 5 times per week.
- Different light modes may be used depending on skin concerns:
 - Red + infrared: anti-aging and radiance
 - Blue: blemishes and oily skin
 - Yellow: soothing and skin tone balancing

CONTRAINDICATIONS

Do not use the LED mask if you:

- Take medications or supplements that may cause light sensitivity (photosensitivity), including certain acne treatments or antibiotics.
- Suffer from photosensitive epilepsy.
- Have a strong sensitivity to light.
- Have eye or retinal conditions.
- Wear a pacemaker or any electronic implant.
- Are pregnant or breastfeeding.
- Have skin lesions, infections, inflammation, active herpes, severe eczema, or psoriasis on the face.
- Have recently undergone intensive cosmetic treatments (laser, peeling, microneedling, injections, etc.).
- Have a known allergy to silicone, plastic, or hook-and-loop fastening materials.
- If you are unsure whether this device is suitable for you, consult your doctor.

PRECAUTIONS FOR USE

- Perform a sensitivity test before first use.
- Start with the lowest intensity setting.
- Clean and dry your skin before each session.
- Do not use on irritated, burned, or injured skin.
- Avoid using photosensitizing skincare products before treatment.
- Do not exceed 30 minutes of use per day.
- Stop using immediately if irritation, discomfort, redness, or a burning sensation occurs.
- Consult a healthcare professional if in doubt.
- If you are undergoing dermatological treatment or have highly sensitive skin, seek medical advice before use.
- Always wear the mask before switching on the LEDs to prevent direct exposure of the eyes to the light.

WARNINGS

- Keep eyes closed during use. Do not stare directly at the illuminated LEDs.
- Risk of Electric Shock. Do not use if the product is damaged.
- For household use only.
- Do not use the mask if the housing, cables, connectors, or battery compartment show signs of damage or deterioration.
- Never disassemble, modify, or attempt to repair the product yourself.
- Do not expose the mask, battery, or charger to water, excessive moisture, or any other liquids.
- Do not use the product while bathing, showering, swimming, or in a highly humid environment.

- Always disconnect the charging cable before using the mask.
- Do not use the mask while it is charging unless explicitly authorized by the manufacturer.
- Use only the charger, cable, and accessories supplied or recommended by the manufacturer.
- Do not insert metal objects into connectors or charging ports.
- If overheating, unusual odor, smoke, or malfunction occurs, discontinue use immediately and disconnect the device from power.
- Do not handle the charger or connectors with wet hands.
- Ensure that the charger power supply matches the product specifications.
- Do not use the product near significant sources of water or liquid splashes.
- Keep the product out of reach of children while charging.
- If the mask falls into water while connected to a power source, do not retrieve it until the power supply has been disconnected.
- Do not use a battery that is swollen, punctured, deformed, leaking, or otherwise damaged.
- Replace any damaged electrical component immediately with a manufacturer-approved part.
- Do not expose the battery to extreme temperatures that could damage its electrical insulation.
- Any maintenance or servicing must be performed by qualified personnel only.
- Failure to follow these instructions may result in electric shock, burns, fire, or serious injury.

Troubleshooting

Issue	Possible Cause	Solution
The mask light is not on or one mode is not working	The USB charging cable is loose	reconnect the USB cable
	Non-matching USB cable	Please use the dedicated USB cable that comes with the mask
	Low battery or depleted	Please charge the controller
Other issues		Please consult the distributor

Note: If the above methods do not resolve the issue, please stop using the device and consult your local distributor.

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