NOOANCE
PARIS

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Table of content

Recommendations for Use	29
Safety Instructions	
Warnings	
Precautions	
Cleaning	
Contraindications	
Potential side effects	
Disclaimer	38
Product Representation	38
Customer Service	38
Product Information	39
Device Specifications	
Contents of the box	
Product Overview	41
Usage	42
Before first use	
First Use and Preparation	
User guide	44
Installation of Sponge Pads	
Frequency of Use	
Storing	
Warranty	
Air travel	
Disposal	
Troubleshooting	48
Symbols	50

Recommendations for Use

Safety Instructions

- Read all instructions carefully before using this device.
- All instructions and this manual are available on the www.nooance-paris.com website or upon request to customer service at: bonjour@nooance.com
- DO NOT charge or use the device near water or if it is at risk of coming into contact with water. Do not immerse the device.
- DO NOT use the device in the bathroom or the shower. DO NOT use if the device has been exposed to water or other liquids. This could lead to an electrical hazard and personal injury.
- DO NOT use the device if it is not working properly, has been dropped, or damaged.
- DO NOT use the device if the timer is faulty. For warranty and repairs, contact customer service.
- DO NOT attempt to repair the device. Repairing or opening the device will void your warranty. The device contains no user-serviceable parts. This device contains non-replaceable batteries.
- DO NOT perform maintenance or disassemble the product. In case of malfunction or damage, contact customer service at bonjour@nooance.com
- DO NOT use the device if it is damaged or the cable is damaged. If any part of the cable is damaged, stop using it. Contact customer service at bonjour@nooance.com to replace the cable.

- DO NOT use or charge the device with a damaged or defective power adapter.
- DO NOT excessively bend or mistreat the USB charging cable. You could damage the device.
- The remote control must only be powered by a safety extra-low voltage of 5V DC.
- Use only the provided USB cable. Using an unsuitable charging cable may damage the device.
- DO NOT wet the helmet or the controller (battery) and cable. This is an electrical hazard that can result in personal injury.
- Please do not use the device while charging the controller(battery). This can lead to faults or electrical hazards.
- DO NOT pull or forcibly detach the power cord when using the device. This could cause a fault or an electrical hazard.
- Keep the device away from hot objects like ovens and indoor heaters.
- If you drop the device or suspect it is faulty(e.g,the LED or laser does not work, the helmet is cracked). Please get in touch with the customer service.

- DO NOT use the device if the cable is damaged, doing so may cause electric shock. Please get in touch with the customer service to replace the defective cable.
- DO NOT attempt to modify or repair the device (helmet, controller).
- Always unplug the device from the power outlet prior to cleaning it. Failure to do so may result in an electrical hazard that may cause harm or injury.
- DO NOT use the device for purposes other than those specified in this manual. Use only the device according to the instruction manual provided.

Warnings

- DO NOT let children manipulate or play with the device. Cleaning and charging the device should not be done by children.
- Before using the device, you must perform a skin sensitivity test.
- DO NOT use the device with wet hair. Make sure your hair is dry before using the helmet.
- Use the device with caution, as handling in a manner other than that specified in the manual may result in a dangerous radiation exposure.
- DO NOT stare directly at the red laser light or its reflection in a mirror because it can temporarily irritate or harm your eyes. Never use a magnifying glass when using the device because it can cause temporary irritation to your eyes.
- DO NOT point the LED light directly into your eyes or another person's eyes. You may experience temporary visual impairment.
- This device can be used by individuals with reduced physical, sensory, or mental abilities, or those lacking experience or knowledge, provided they are properly supervised or have been instructed on the safe use of the device and understand the risks involved.
- STOP using the device if it overheats or if your skin emits an uncomfortable warmth. Turn off the device immediately and contact the customer service.

- DO NOT use abrasive solvents or abrasive cleaners other than 70% alcohol on the device.
- DO NOT use the device excessively – Respect the recommended conditions of use and frequency of use.
- DO NOT use the device with unauthorized accessories.
- No modifications to this equipment are allowed.
- DO NOT use the device for purposes other than those specified in this manual. Use only the device according to the instruction manual provided.
- Contact your doctor if you have any medical questions regarding the use of this device.

Precautions

- DO NOT scratch or damage the light-emitting surface.
- Your warranty will be void if you do not use the device in accordance with the manual.
- DO NOT use the device outdoors to avoid damage.
- DO NOT place the device in direct sunlight.
- DO NOT sit, lie down, or step on the device to avoid breaking it.
- DO NOT use or store the device near hot surfaces.
- Store the device according to the protocol described in this manual.
- The device is designed for cosmetic use only by healthy individuals. Please contact your doctor if you are unsure of your health condition.
- The device's battery is not removable by the user. If the device needs to be disposed of, do not burn it. Dispose of the device in accordance with local regulations.
- Inspect the device and do not use it if you see signs of damage, such as battery leaks. Contact our customer service.
- We strongly recommend exercising caution when taking medications. If in doubt, consult your doctor.

Note: For hygiene reasons, we advise not sharing your device with others without disinfecting the part in contact with the scalp skin using 70% alcohol. Never use acetone or other solvents on any part of your device. This can block the laser light from coming through.

Cleaning

- Clean the device with a soft, dry cloth.
- DO NOT use abrasive solvents or abrasive cleaners other than 70% alcohol on the device.
- The device should not be cleaned or charged by children.
- The helmet is reusable, and it needs to be cleaned after each use. Follow the instructions below to clean your device.
 - 1)Before cleaning, please turn off the device and unplug the power.
 - 2)Use a clean cotton cloth moistened with water to clean the casing and ensure that no liquid enters the interior of the device. Stubborn stains can be removed by wiping gently with 70% alcohol.
 - 3)Visually check the cleanliness of the device, make sure there is no dust, impurities, or any other pollutant, otherwise, repeat the above step again.
 - 4)Drying: dry the device in the air naturally.

Note: Never immerse the device in water or other cleaning liquids to avoid damage to the device.

Contraindications

Do not use the device if you are concerned by any of the following situation:

- DO NOT use the device if you are under the age of 18.
- DO NOT use the device if the user is pregnant or think you may be, or breastfeeding.
- DO NOT use the device if you have precancerous lesions, cancerous lesions or active cancer.
- DO NOT use the device if you are photosensitive.
- DO NOT use the device is you are taking photosensitive drugs.
- DO NOT use the device if you have a medical history of seizures that were triggered by light.
- DO NOT use the device if you have skin problems such as lupus erythematosus, photosensitive eczema, albinism, or cancerous lesions as you may cause a severe skin reaction.
- DO NOT use if you have light induced headaches or migraines.
- DO NOT use the device if you have head's open wounds
- DO NOT use the device if you have had recent surgery or other surgical procedures.
- DO NOT use the device on bacterial infections as this could worsen the symptoms.

Potential side effects

- Normal Reactions: You may feel a slight tingling or warm sensation during or after the session. Your skin may appear slightly red. These effects should subside within 1 to 2 hours.
- When to Stop the use of the device / when to contact the customer service (bonjour@nooance.com) and/or consult a Doctor :
 - If redness, itching, or tingling persists more than 2 hours after a session.
 - If you develop a persistent headache or migraine.
 - If the device feels uncomfortably hot or causes a burning sensation (Stop immediately).
 - If you experience any other unexpected reaction within 48 hours of your first session.

Disclaimer

- This device is not intended to diagnose, treat, or prevent any disease. NOOANCE makes no claims or warranties regarding this product's ability to treat any mental, physical, or dermatological condition. Consult a qualified healthcare professional for any condition requiring medical attention.
- The device is not intended for use in the diagnosis of diseases and other conditions, nor for the cure, mitigation, treatment, or prevention of diseases, and is not intended to alter the structure or function of the body.
- Aesthetic, non-medical device intended for home use.

Product Representation

NOOANCE continually improves its products using the latest technologies. As a result, the actual product may differ slightly from the product shown in the images on the packaging and in this user guide.

Customer Service






For more information on using or maintaining the device, or to obtain additional information, please contact Customer Service via email at bonjour@nooance.com.

Product Information

Device Specifications


Product Name	Laser and LED HELMET X400
Model	M-282A, M-282B
LED output optical power density	<5mW / cm ²
Laser power density	<5mW
LED and laser wavelengths	655 nm±5nm
Number of Light Emitting Diodes(LED)	300
Number of Lasers	100
Laser classification	Class 1
Mode	Single-mode continuous mode with a 22Khz pulse
Operating temperature	0°C to 40°C 32°F to 104°F
Storage and transport temperature	- 5°C to 55°C 23°F to 131°F
Operating humidity	30% to 95% without condensation
Storage and transport humidity	20% to 95% without condensation
Charging supply	5V $\overline{=}$ 2A
Charging port	USB-C
Charging cable	USB-C to USB-C
Charging cable length	1 meter / 3.28 feet
Battery	3.7V, 5200mAh
Battery type	Lithium-ion
Charging time	240 minutes/ 4 hours
Recommended use	3 times per week
Timer	15 minutes
Dimension of the device	243.38*213.88*172.14mm
Weight of the device	0.652kg / 23oz

Contents of the box

No.	Accessories	Quantity
1	Laser and LED HELMET X400 (M-282A, M-282B)	1
2	Sponge Pad (Size: 0.5cm, 1.5cm)	6
3	User manual	1
4	Controller (battery)	1
5	USB Cable	1
6	Storage Bag	1

Product Overview

- ①. Display screen
- ②. Power ON/OFF button
- ③. Pause / Start button
- ④. Charging port
- ⑤. USB-C port
- ⑥. Laser emission
- ⑦. LED emission

Usage

Before first use

We strongly recommend that users carry out a light sensitivity test before first use.

1. Wear the device properly on your head, ensuring it is correctly positioned and in full contact with the scalp. The device will only activate when worn.

2. Switch on the device to start the session and allow it to operate for 5 minutes.

3. Turn off the device after 5 minutes (the default timer is set to 15 minutes).

4. After use, remove the device and visually inspect the scalp area.

Visually inspect the treated area. If your skin feels excessively warm or itchy during the session, stop immediately. If the area turns red and this redness persists for more than 2 hours, your skin may be hypersensitive to light. Discontinue use and contact customer service at: bonjour@nooance.com.

First Use and Preparation

1)Unpack the package and check off all the parts described in section « Contents of the box »

2)Check the charge of the device / How to charge the device:


- Your device may be partially or fully charged upon arrival. To verify the battery level, insert the provided USB-C cable on the remote control (you should hear two clicks), then connect the cable to a safety extra low voltage 5V adaptor.
- A steady indicator light confirms a sufficient charge: you can use your device without further charging it.
- A blinking light indicates that the device requires further charging. Keep the remote control in charge.
- Once charging is complete (after 4 hours), the light will remain steady.
- The battery should be recharged on average every 8 to 10 sessions.

When the controller battery level is less than 10%, the indicator lights flash white. To recharge the device, follow the instructions above.

User guide


Preparation: Make sure your hair is dry before using the helmet.

Step1

Please select the appropriate sponge pad size to fit inside the helmet.


Step 2

Connect the helmet to the controller (battery) using the provided cable : you should hear 2 clicks

Step 3


Position the helmet in its upright position on the top of your head.

Step 4

Press and hold the ON/OFF button for 3 seconds to switch the device on.


(Note: If the start/pause button is not pressed within five minutes, the device will turn off automatically.)

Step 5

Press quickly the "start/pause" button (WITHOUT holding it) to begin the 15-minute session, the light turn on indicates the session has started. Once the session is complete, the helmet light will turn off automatically.

IMPORTANT: Unless the NOOANCE helmet is positioned on the user's head, it will not begin a session. If the NOOANCE helmet is removed from the head during a session, it will pause automatically. To resume the session, position the NOOANCE helmet on the user's head within 5 minutes and press the start/pause button.

Step 6

If you need to pause the session at any time, press the Start/Pause button and the timer will stop.

Press the Start/Pause button again to resume the session. If the start/pause button is not pressed within 5 minutes during a paused session, the NOOANCE helmet will turn off automatically.

Step 7

To turn off the helmet, hold the button for 3 seconds.

IMPORTANT: The controller(battery) has a memory in case of a power failure. In the event of a power failure or accidental disconnection of the USB cable during use, simply press the start/pause button again. The helmet will resume the timer from where it left off. To switch off the device, press the ON/OFF button for 3 seconds. Then follow steps 2 to 7.

Installation of Sponge Pads



There are 2 sizes to choose from. The corresponding number is located below each pad. Please put the correct number in the specified position.

Please note: contact between the helmet and the scalp is normal. The sponge pads are designed to stabilise the helmet during use.

0.5CM

1.5CM

Frequency of Use

We recommend using the helmet for 15 minutes, 3 times a week, during 16 weeks.

Storing

- 1) Avoid direct sunlight on the device while storing.
- 2) It is recommended to store the device in its original packaging.
- 3) Store in an environment with a temperature of -5°C / $+55^{\circ}\text{C}$, a humidity of 20-95%RH and an atmospheric pressure of 700hPa-1060hPa.

Warranty

- The device comes with a 24-month guarantee for the original purchaser, subject to proof of purchase from an authorised vendor. This guarantee covers normal use with approved accessories. In the event of a fault, contact our customer service department.
- Please note that the guarantee does not cover damage caused by accidents, incorrect use, use of non-recommended accessories, or poor electrical connection.

Air travel

The device is allowed in the cabin. If you wish to check it into the hold, take the controller/Battery with you in the cabin. The main airlines refuse to accept lithium batteries in the hold. If in doubt, consult your airline.

Disposal

DO NOT dispose of the device with household waste. The device should be recycled. The device contains a Lithium-ion battery. For protection of the environment, please safely dispose of the device following your local guidelines.


Troubleshooting


If the actions suggested fail to fix the problem, or the problem reoccurs, and your device is within the warranty period please contact our customer service by writing to bonjour@nooance.com.

- DO NOT attempt to investigate any further.
- DO NOT attempt to modify the device in any way. The device has NO parts you can repair yourself.

Problem	Cause	Solution
The device doesn't switch on	The controller(battery) is not fully charged because the charging cable was not properly inserted into the power outlet	Insert the charging cable properly into the power outlet and device
	The user does not wear the helmet upright on the head.	Unless the helmet is positioned on the user's head , the session will not begin.
	Allow the helmet internal software to reset itself	Hold for 3–5 seconds and allow the helmet's internal software to reset itself.
The helmet is powered on but pressing does not activate the helmet.	The helmet will not power on unless the electronic sensor detects that the helmet is positioned on the user's head. This is a built-in safety feature to protect the eyes from direct laser light exposure.	Activate the helmet by placing it vertically over your head.
The helmet is not Lighting up, but the controller is on	The power of the controller is low.	Fuly charge the controller before use.
	If the device still fails to light up after following the instructions and troubleshooting.	Please get in touch with customer service
The helmet controller doesn't show the numbers clearly on the timer.	Try shutting down the device, restarting it, and allowing the helmet's internal software to reset itself.	If the timer cannot be displayed clearly after the helmet is restarted, please get in touch with customer service.

Symbols

symbol	Explanation
	Indicates that the device is CE marked and authorized for sale in Europe in accordance with applicable safety and conformity assessment standards (refer to the EU Declaration of Conformity, available upon request, for the list of relevant standards).
	The Waste of Electrical and Electronic Equipment (WEEE), indicates the device may not be discarded as household waste (2012/19/EU).
	Indicates that the electromagnetic radiation from the device is below the limits specified by the Federal Communications Commission and the manufacturer has followed the requirements of the Supplier's Declaration of Conformity authorization procedures.
	This warning symbol on the device label is a warning to read the enclosed documents before use and take all safety precautions required.
	Consult instructions for use.
	Date of Manufacture.
	This symbol means keep the device dry.
	This symbol indicates the serial number of the device.
	This symbol refers to the batch code of the device.
	This symbol indicates the USB connector.
	This symbol means that the device uses a Lithium-ion Polymer battery.
	This symbol means keep the device away from direct sunlight.
	This symbol means the upper and lower operating temperatures that the device can be safely used at.

	This symbol indicates the upper and lower operating humidity that the device can be safely used at.
	This symbol means that the cardboard packaging that the device came in can be recycled.
	The BF symbol, which indicates this device is in accordance to the degree of protecting against electric shock for type BF equipment.
	Protected against solid foreign objects of 12.5 mm and greater, and Protection against vertically falling water drops when ENCLOSURE tilted up to 15°.
	Unique Device Identifier.
	Indicates that the device is UKCA marked and authorized for sale in the United Kingdom in accordance with applicable safety and conformity assessment standards (refer to the UKCA Declaration of Conformity, available upon request, for the list of relevant standards).
	Indicates that the device is tested according to RoHS requirements.
	Warning: laser radiation
	Manufacturer
	Country of manufacture
	Importer
	EU Representative