

PHILIPS sonicare

For Kids 3+

English 6
Azərbaycanca 21
Čeština 37
Dansk 52
Deutsch 67
Eesti 86
Español 102
Français 119
Hrvatski 136
Italiano 152
Latviešu 169
Lietuviškai 184
Magyar 200
Nederlands 217
Norsk 234
O'zbek 249
Polski 265
Português 282
Română 299
Shqip 315
Slovenščina 331
Slovensky 346
Srpski 362
Suomi 378
Svenska 394
Türkçe 409
Türkmençe 425
Ελληνικά 441
Български 459
Кыргызча 477
Македонски 494
Русский 511
Тоҷикӣ 530

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Introduction

Congratulations on your purchase and welcome to Philips! To fully benefit from the support that Philips offers, register your product at www.philips.com/welcome.

Important safety information

Only use the product for its intended purpose. Read this information carefully before you use the product and its batteries and accessories, and save it for future reference. Misuse can lead to hazards or serious injuries. Accessories supplied may vary for different products. NOT A TOY. Check brush heads regularly, cracks from dropping can lead to a choking hazard.

WARNINGS

- Only use original Philips accessories or consumables. If supplied, only use the original USB cable and USB power supply unit.
- This appliance contains no user-serviceable parts. If the appliance is damaged, stop using it and contact the Consumer Care Center in your country (see 'Warranty and support').
- Keep charging cables and wall adapters away from water.
- After cleaning, make sure the wall adapters and charging cables are completely dry before connecting to power.
- This appliance can be used by children and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Cleaning and user maintenance shall not be made by children without supervision.
- Children should be supervised to ensure that they do not play with the product.
- Do not clean any part of the product in the dishwasher unless otherwise instructed in the Cleaning section.

- Do not use charging cables and wall adapters outdoors or near heated surfaces.
- If the charging cable or wall adapter is damaged, stop using it.
- This appliance has only been designed for cleaning teeth and gums.
- Replace the brush head every 3 months or sooner if signs of wear appear. Stop using a brush head with crushed or bent bristles.
- Rinse the handle, especially the brush head connection. Gently clean around the rubber seal.
- Make sure brush head and toothbrush are dry before storing in travel case.
- Avoid direct contact with products that contain essential oils or coconut oil. Contact may result in bristles dislodging.
- The maximum altitude of use is 4500 meters.
- Information on REACH substances, contained in Philips products, can be found on www.philips.com/REACH

MEDICAL WARNINGS

- Stop using the appliance and consult your dentist/doctor if excessive or prolonged bleeding occurs after use, or if you experience discomfort or pain during use.
- If you have a pacemaker or other implanted device, contact your doctor or the manufacturer of the implanted device prior to use.
- This appliance is a personal care device and is not intended for use on multiple patients in a dental practice or institution.

Battery safety instructions

- This appliance contains non-replaceable batteries that are only to be removed by skilled persons.
- Keep product and batteries away from fire and do not expose them to direct sunlight or high temperatures.
- If the product becomes abnormally hot or smelly, changes color or if charging takes longer than usual, stop using and charging the product and contact Philips.

- Do not place products and their batteries in microwave ovens or on induction cookers.
- Do not open, modify, pierce, damage or dismantle the product or battery to prevent batteries from heating up or releasing toxic or hazardous substances. Do not short-circuit, overcharge or reverse charge batteries.
- This appliance contains batteries that are non-replaceable. When the battery is at end of life, the appliance shall be properly disposed of, see section on Recycling.
- Charge only with Philips WAA1001, WAA2001 safety approved wall adapters with input ratings (100 -240V~; 50/60Hz; 3.5W) and output ratings (5V ; 2.5W). For the list of approved adapters, please see **www.philips.com/support**.
- Charge and use the product at a temperature between 0 °C and 40 °C (32 °F and 104 °F).

Storage and transportation conditions

Temperature -20 °C to +60 °C (-4 °F to 140 °F).

Electromagnetic fields (EMF)

This Philips appliance complies with all applicable standards and regulations regarding exposure to electromagnetic fields.


LIST OF SYMBOLS USED

The following symbols may appear on the product:

This symbol means read the instructions for use, before using the appliance.

This symbol means to use the indicated power supply.

Intended Use

Sonicare for Kids power toothbrushes are intended to remove adherent plaque and food debris from the teeth to reduce tooth decay and improve and maintain oral health. Sonicare for Kids power toothbrushes are intended for personal use. Use by children should be with adult supervision.

Your Philips Sonicare toothbrush

- 1 Handle
- 2 Power on/off button
- 3 Mode indicator
- 4 Battery indicator
- 5 USB charger (wall adapter not included)

Accessories:

- 6 Brush head cap
- 7 Brush head
- 8 Stickers
- 9 Travel case

Note: Accessories included may vary based on the model purchased.

Getting started

Sonicare for Kids App

Download the Sonicare for Kids App to support healthy habit building. The App is compatible with wide range of iOS and Android^(TM) smart devices.

Note: If your toothbrush is not Bluetooth[®] capable you can still use the Sonicare for Kids app on iOS and Android^(TM) smart devices to enhance your brushing experience.

Applying your favorite sticker

This toothbrush comes with several stickers to choose from.

- 1 Apply sticker(s) on the handle to create your personalized toothbrush.

Note: Apply stickers only on the handle and travel case. Do not apply stickers on the brush head.


Note: Make sure the handle or travel case is clean and dry, and that stickers are fully pressed onto the handle with no air gaps to ensure full adhesion.

Note: Kids should not be left unattended when applying the stickers.

Attaching the brush head

- 1 Align the brush head so the bristles point in the same direction as the front of the handle.

- 2 Firmly press the brush head down onto the metal shaft until it will not go any further.

Note: It is normal to see a slight gap between the brush head and the handle.


Using your Sonicare for Kids toothbrush

If you are using your Sonicare for Kids toothbrush for the first time, it is normal to feel more vibration than when using a non-electric toothbrush. It is common for first-time users to apply too much pressure. Only apply gentle pressure and let the toothbrush do the brushing for you. Please follow the brushing instructions below for the best experience.


Your Sonicare for Kids toothbrush comes with the EasyStart feature activated to help you transition to using a power toothbrush. This feature gently increases the power over the first 14 brushing sessions to help you adapt to brushing with a Philips Sonicare toothbrush.

Note: For children under age 8, it is recommended that a parent or adult caregiver perform or supervise toothbrushing.

Brushing instructions


- 1 Wet the bristles and apply a small amount of toothpaste.


- 2 Place the toothbrush bristles against the teeth at a slight angle (45 degrees), pressing gently to make the bristles reach the gumline or slightly beneath the gumline.

Note: Keep the center of the brush in contact with the teeth at all times to ensure the whole tooth surface is brushed.

- 3 Press the power on/off button to switch on the toothbrush.

- 4 Gently keep the bristles placed on the teeth and in the gumline. Brush your teeth with small back and forth motion so the bristles reach between the teeth. Continue this motion throughout your brushing cycle.

Note: The bristles should slightly flare. It is not recommended to scrub as you would with a manual toothbrush.

- 5 To clean the inside surfaces of the front teeth, tilt the brush handle semi-upright and make several vertical overlapping brushing strokes on each tooth.

Note: To make sure you brush evenly throughout the mouth, divide your mouth into 4 sections using the Sound Pacer feature (see chapter 'Features').


- 6 After you have completed the brushing cycle, you can spend additional time brushing the chewing surfaces of your teeth and areas where staining occurs.

Note: Your Philips Sonicare toothbrush is safe to use on braces (brush heads wear out sooner when used on braces) and dental restorations (fillings, crowns, veneers) if they are properly adhered and not compromised. If you experience any issues, consult your dental professional.

***Disclaimer**

If the Philips Sonicare for Kids toothbrush is used in clinical studies, it must be used in Gentle mode with the handle fully charged and the EasyStart feature deactivated.

Modes

Your power toothbrush is equipped with two different modes: Extra Gentle, and Gentle.


To change the mode:

- Press the power button once to turn on the toothbrush. Press a second time within **2 seconds** to change the brushing mode.
- After 2 seconds of brushing, you will need to press the power button to pause first, and then press it again **within 2 seconds** to switch modes.


Mode	Extra Gentle	Gentle
Benefit	Reduced power mode used for training and introducing younger children to the sonic experience.	Full-strength mode designed for older/independent brushers.
Total brushing time	2 minutes	2 minutes

How to brush

Brush each segment for 30 seconds

Brush each segment for 30 seconds

Features

- 2-minute Sound Timer
- Sound Pacer
- EasyStart

2-minute Sound Timer

The 2-minute Sound Timer plays a celebratory sound at the end of the predetermined brushing time. Additionally, the timer automatically turns off the handle when the brushing cycle is complete.


- In Extra-Gentle Mode, the timer is set at 1 minute and slowly increases over time. This gradual increase helps young children build up to the dentist-recommended 2 minutes of brushing.
- In Gentle Mode, the timer is set at 2 minutes.

Dental professionals recommend brushing no less than 2 minutes twice a day.

Note: Instruct your child to brush until they hear the celebratory tune and the toothbrush turns off.

Note: If you press the power on/off button after you started the brushing cycle, the toothbrush pauses. If you press the on/off button within 30 seconds, the timer resumes to allow you to finish your brushing cycle. After a pause of 30 seconds, the timer resets.

Sound Pacer

To make sure your child brushes evenly throughout the mouth, divide the mouth into 4 sections using the Sound Pacer feature.

The pacer splits brushing time into 4 equal segments and signals with a short tune when it's time to move to the next area.

Note: Instruct your child to brush until they hear the celebratory tune and the toothbrush turns off.


EasyStart

This Philips Sonicare model comes with the EasyStart feature activated. The EasyStart feature gently increases the power over the first 14 brushing sessions to help you get used to brushing with the Philips Sonicare.


Note: Each of the first 14 brushing sessions must be at least 1 minute in length to properly advance through the EasyStart ramp-up cycle.

Activating or deactivating the EasyStart feature

Step 1: Place the handle on a plugged in charger.


Step 2: Press and hold the power on/off button until you hear a single short beep (for 3 seconds).

If you see the battery indicator flash green 2 times and hear 3 tones low to high, then the feature has been activated.

OR

If you see the battery indicator flash amber 1 time and hear 3 tones high to low, then the feature has been deactivated.


Step 3: Release the power on/off button.

Charging and battery status

This Philips Sonicare toothbrush is designed to provide on average 42 brushing sessions, each session lasting 2 minutes (21 days if used twice per day) on a fully charged battery. The battery indicator shows the battery status when you complete the 2-minute brushing session, when you pause the toothbrush, or while it is charging.

Note: It can take up to 24 hours to fully charge the battery before first use, but you can use the Philips Sonicare toothbrush before it is fully charged.

Charging

- 1 Attach the USB charging cord to a wall adapter, plug the wall adapter into an electrical wall outlet.
- 2 Place the handle on the charger.
- 3 Two short beeps are emitted to confirm that the handle is placed properly.
- 4 The flashing green light of the battery indicator shows that the toothbrush is charging.
- 5 Once the battery is fully charged, the green light remains solid for 30 seconds and then turns off.

Battery status (when handle is on plugged-in charger)

When the handle is plugged into the charger, the battery indicator at the bottom of the handle shows the battery level.

Battery status	LED color
Charging	Flashing green

Full

Solid green for 30 seconds and then switches off

Battery status (when handle is not on plugged-in charger)


When the toothbrush is awake, the battery indicator at the bottom of the handle shows the battery level.

Battery status	LED color
Full	Solid green
Medium low	Flashing green
Low	Flashing amber and 3 beeps after brushing
Empty	Flashing amber and 2 sets of 5 beeps

Note: To save energy, the battery indicator will switch off when not in use.


Note: When the battery is completely empty, the toothbrush switches off. Place the toothbrush on the charger to charge it.

Note: To keep the battery fully charged at all times, you may keep your toothbrush plugged into the charger when not in use.

Cleaning**Toothbrush handle**

- 1 Remove the brush head and rinse the metal shaft area with warm water.

Caution: Do not push on the rubber seal on the metal shaft or around the buttons with sharp objects, as this may cause damage.


- 2 Make sure you remove any residual toothpaste on the handle and around the buttons at least once a week, but preferably after each use.
- 3 Wipe the entire surface of the handle with a damp cloth.

Caution: Do not tap the handle on the sink or counter to remove excess water.

- 4 Dry the handle with a soft cloth or towel.

Brush head

- 1 Remove the brush head from the handle and rinse the brush head connection with warm water at least once a week.
- 2 Rinse the brush head and bristles after each use and allow to air dry.

Charger

- 1 Unplug the charger before you clean it.
- 2 Wipe the surface of the charger with a damp cloth.
- 3 Dry the charger with a soft cloth or towel.

Travel case

Rinse with warm water and use a damp cloth to clean the travel case and allow to air dry.

Storage

If you are not going to use the product for an extended period of time, unplug it from the electrical outlet, clean it and store it in a cool and dry place away from direct sunlight.

Replacement

Replace Philips Sonicare brush heads every 3 months to achieve optimal results. Use only Philips Sonicare replacement brush heads.

Locating the model number

The model number begins with "HX" (e.g. HX123A) and can be found on the bottom of your toothbrush.

Troubleshooting

This chapter summarizes the most common problems you could encounter with the appliance. If you are unable to solve the problem with the information below, visit **www.philips.com/support** and insert your model number for a list of frequently asked questions or contact the Consumer Care Center in your country.

Problem	Possible Cause	Solution
I cannot fully attach the brush head. There is a gap between the brush head and handle.		This gap is normal and is needed for the brush head to vibrate appropriately. The brush head must be able to move to create the right amount of vibrations.
My Philips Sonicare toothbrush vibration is less powerful than before.	The brush head is too close to the body of the toothbrush or the toothbrush needs to be charged.	Remove the brush head from the handle and re-set it on the metal shaft leaving a small gap. (1-2 mm). When charging, make sure the charger is plugged into a power supply and placed on a flat surface before placing the toothbrush on the charger. Allow the toothbrush to charge for 24 hours.

Problem	Possible Cause	Solution
I charged my toothbrush, but it only works for a short time.	The base of the toothbrush is not placed flat on the charger when attempting to charge.	Make sure the charger is plugged into a power supply and placed on a flat surface before placing the toothbrush on the charger. Allow the toothbrush to charge for 24 hours.

Warranty and support

Please visit **www.philips.com/support** and enter the model number located at the bottom of your device in the search bar in order to obtain more information on how to register your device, warranty information, and other frequently visited support topics.

For all product support, power consumption details, and the EU Declaration of Conformity, please visit **www.philips.com**. After searching for your device, you can find information on your product's features, technical specifications, as well as available accessories.

Warranty exclusions

What is not covered under warranty:

- In mouth attachments including brush heads and nozzles.
- Damage caused by use of unauthorized replacement parts or unauthorized brush heads.
- Damage caused by misuse, abuse, neglect, alterations or unauthorized repair.
- Normal wear and tear, including chips, scratches, abrasions, discoloration or fading.

Recycling

- This symbol means that electrical products and batteries shall not be disposed of with normal household waste.
- Follow your country's rules for the separate collection of electrical products and batteries.

Removal of built-in rechargeable battery

- CONTAINS LITHIUM-ION RECHARGEABLE CELL. MUST BE RECYCLED OR DISPOSED OF PROPERLY.
- The built-in rechargeable battery can be removed for disposal. If you are not comfortable disassembling electronics, please seek professional assistance at your local electronics recycling or repair shop.
- Instructions for removal of built-in rechargeable batteries can be found on **www.philips.com/support**. Use the Philips support home search bar to enter the toothbrush model number and find your product. The model number begins with "HX" (e.g. HX123A) and can be found on the bottom of your toothbrush. The battery removal instructions are in the Troubleshooting and Repair section of the product's support page.
- When handling batteries, make sure that hands, the product and the batteries are dry.
- If batteries are damaged or leaking, avoid contact with the skin or eyes. If this occurs, immediately rinse well with water and seek medical care.
- To avoid accidental short-circuiting of batteries after removal, do not let battery terminals come into contact with metal objects (e.g. coins, hairpins, rings). Do not wrap batteries in aluminum foil. Tape battery terminals or put batteries in a plastic bag before you discard them.

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© 2024 Koninklijke Philips N.V. (KPNV). All rights reserved. Philips and the Philips shield are trademarks of KPNV. Sonicare and the Sonicare logo are trademarks of Philips Oral Healthcare, LLC and/or KPNV.
The Bluetooth®/MD word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Koninklijke Philips N.V. is under license.