medisana®


LM 200

medisana GmbH
Carl-Schurz-Str. 2,
41460 NEUSS,
Germany

REF **88401**

EAN **40 15588 88401 6**

Kompressions-Bein-Massagegerät LM 200 (DE)

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Air Compression Leg Massager LM 200 (GB)

❶ Massager with control panel
❷ Indicator for battery level ❸ On/Off button
❹ Level-button, with indicator above:
- » low intensity – the indicator is lit in green colour.
- » medium intensity – the indicator is lit in blue colour.
- » high intensity – the indicator is lit in red colour.

❺ Heat-button, with indicator above:
- » low heat level – the indicator is lit in green colour.
- » medium heat level – the indicator is lit in blue colour.
- » high heat level – the indicator is lit in red colour.

❻ Number display for the remaining operating time, and for the set intensity and the set heat level
❼ Charging port (USB connector)
❽…⓫ How to fasten the massager onto the lower leg.

Safety Instructions

Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use. These brief instructions are an integral part of the instructions for use. The complete instructions for use are available for download at https://www.medisana.com. Should you give this device to another person, it is vital that you also pass on the instructions for use.

- Please ensure that the polythene packing is kept away from the reach of children! Risk of suffocation!
- The device must only be used for its intended purpose as described in the instruction manual. Using the device for any other purpose invalidates the warranty.
- The device is not to be used for commercial use or medical purposes. If you have health concerns, consult your doctor before using the device.
- Only use the massager indoors.
- Do not use the massager in damp rooms (e.g. when bathing or showering).
- Make sure that the maximum operating time of 15 minutes is not exceeded!
- Keep the device away from direct sunlight, do not expose it to high temperatures.
- Follow all legal constraints and regulations for disposing of of the device and the internal battery. (See also Section “Disposal”.)
- This device must only be supplied at safety extra low voltage (SELV) corresponding to the marking on the appliance.
- Do not carry, hold or drag the device by the cable!
- Batteries must never be squeezed, disassembled or affected by hot temperatures (above 60°C). There is the danger of fire, explosions and the risk of a leakage of dangerous and poisonous chemicals.
- This device contains non-replaceable batteries.
- The batteries contained in this device may be replaced only by service technicians that have been authorised by **medisana**.
- The battery life gets reduced if it is not used for a long time. It is therefore recommended to charge the battery at least once or twice per month.
- Use a USB mains adapter with an output of 5V and (at least) 1A for charging the battery. A larger output voltage will damage the device.
- Only insert the power supply unit into a socket when the device is switched off. The device must not be used while it is being charged.
- Do not touch the mains adapter if you are standing in water or if your hands are wet or damp. Otherwise, there might be the risk of damages or an electric shock.
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children must be supervised to ensure that they do not play with the device.
- Do not use this unit to supplement or replace medical treatment. Chronic complaints and symptoms may be made worse.
- Do not use the device, or only after consulting your doctor, if:
 - » you are pregnant,
 - » you have a cardiac pacemaker, artificial joints or electronic implants,
 - » you suffer from one or more of the following illnesses and complaints: Circulation disorders, varicose veins, open wounds, bruises, broken skin, inflammation of the veins.
- If you experience any pain or discomfort during the massage, stop immediately and consult your doctor.
- The surface of the appliance is hot. People who are insensitive to heat must be careful when using the appliance.
- Before using the massager, consult your doctor regarding any unexplained pain, if you are undergoing medical treatment or using medical equipment.
- Children under the age of three must not use this device due to their inability to react to overheating.
- Always carefully inspect the device for damage before use. Do not attempt to use a defective device.
- Do not use the device if it shows signs of damage, if it is not in perfect working order, or if it has been dropped or become damp. In order to prevent hazards, always send the unit in to the service centre.
- Only use the device in the proper place, that is, around the neck and on the shoulders.
- Do not use the device while you are asleep or when you are lying in bed.
- Do not use the device before going to sleep. The massage has a stimulating effect.
- Never use the device while driving a vehicle or operating a machine.
- The maximum operating time for one session is 15 minutes. Prolonged use reduces the service life of the device.
- Very long periods of use can lead to overheating. Always allow the device to cool down before using it again.
- Never put the unit down next to an electric stove or other sources of heat.
- Do not allow the massage cushion to come into contact with pointed or sharp objects.
- Do not cover up the appliance when it is switched on. Never use it under blankets or cushions. There is a risk of fire, electric shock or injury.
- The device is maintenance free. Do not try to open it. There are no parts inside that are to be cleaned or serviced by the user.
- You may only clean the device yourself. If the device malfunctions please do not repair it yourself as this leads to the automatic expiry of all warranty rights. Consult your specialist dealer and have repairs carried out by the service centre to avoid dangers.
- Do not immerse the device in water or other fluids. If, nonetheless, liquid should penetrate the device, unplug the mains cable immediately.
- Cleaning and user maintenance must not be carried out by children without supervision.

Charging the battery

- Use a 5V mains adapter with an output current of (at least) 1 A. Connect the supplied USB charging cable to the adapter and to the USB port ❼ of the device.
- **The device cannot be used during the charging process.**
- The charging process is indicated by the flashing of one bar of the battery level indicator ❷. When the battery is fully charged, alls bars of the battery level indicator ❷ are lit continuously (without flashing).
- If the lowest bar of the battery level indicator ❷ is flashing the remaining charge is low – the battery should be recharged as soon as possible.

Using the device

Fasten the unit onto your lower leg, as shown in illustrations ❽ to ⓫. Make sure that it is fastened safely and confortably, and that the orientation is correct – the control panel must face outwards.

- Switch the unit on with the ON/OFF button ❸. The massage starts with low intensity and without heat. The number display ❻ is counting downwards the remaining minutes. The maximum operating time is 15 minutes, afterwards the device switches off automatically.
- You can switch between three different intensities using the level button ❹. The currently active intensity (“P1”/“P2”/“P3”) gets displayed after each button press for a few seconds. Moreover, the intensity is visible by the colour of the indicator that is located above the button ❹: green → blue → red.
- Using the heat button ❺, the heating function can be switched on in 3 stages. The currently active heat level (“H1”/“H2”/“H3”) gets displayed after each button press for a few seconds. Moreover, the heat level is visible by the colour of the indicator that is located above the button ❺: green → blue → red.
- You can switch off the appliance at any time by pressing the ON/OFF button ❸. After 15 minutes, the appliance switches off automatically. Then let the appliance cool down for at least 15 minutes before using it again.
- The inner cover is fixed via four touch fasteners, so that it can be simply detached and washed gently by washing machine at 30°C. Afterwards, let it dry by hanging on a line. It is not permitted to use a tumble drier!
- The massager (i.e. the device itself) may be cleaned only with a slightly damp sponge. Never use brushes, aggressive cleaning agents, petrol, thinners or alcohol. Rub the unit dry with a clean dry cloth.

Voltage supply: Li-Ion battery 3.7V, 2500mAh,
Input: 5V⎓ 1.0A

Disposal: This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Consult your local authority or your supplier for information about disposal.

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