

# Braun ThermoScan® Hygiene caps

LF40EULA01



For accurate and more hygienic temperature monitoring.

Braun disposable hygiene caps enable accurate and hygienic temperature taking. It is important to replace hygiene caps after each measurement. The disposable hygiene caps ensure hygienic use without having to clean or disinfect the sensor and virtually eliminate cross contamination. What's more, fitting a new Braun hygiene caps each time ensures that dirt or ear wax do not build up and clog the sensor window, which could lead to inaccurate temperature readings.

Hygiene caps can be used with all models of Braun ThermoScan® IRT ear thermometers.

**Includes**

40x Hygiene caps

## Key Product Highlights



### Ear Thermometer Compatibility

Suitable for all Braun ThermoScan® ear thermometers.



### Safe & Hygienic

Hygienic caps protect the thermometer lens from scratches and impurities which might affect the readings.



### No Cross-Contamination

Eliminates the risk of transferring germs by cross-contamination.



### Not Made of BPA & Latex


The product does not contain BPA plastic or Latex.



## Technical information LF40EULA01

	Quantity	Weight (kg)	Dimensions (mm)	EAN Code	PZN
<b>Unit box</b>		0.05	52 x 34 x 12	04022167400062	7437651
<b>Inner case</b>	6	0.23		24022167400066	
<b>Master case</b>	432	1.689	64 x 79 x 148	14022167400069	
<b>Nb. layers</b>	EU : 4 UK: 2				
<b>Nb. cases per layers</b>	EU : 2 UK: 3				
<b>Pallet</b>	EU: 3456 UK: 2592		EU: 120 x 80 x 195 UK: 57.5 x 74 x 44.5		

## Braun Healthcare



ThermoScan® is a registered trademark of Helen of Troy Limited and/or its affiliates.  
 Certain trademarks used under license from The Procter & Gamble Company or its affiliates.  
 © 2023, All rights reserved.



**Kaz Europe Sàrl**  
**A Helen of Troy Company**  
 Q-Center, Chaux 4  
 CH-1030 Bussigny - Switzerland  
 Phone: +41 21 644 0110  
 Fax: +41 21 644 0111

Website: [www.HelenofTroy.com/emea-en/](http://www.HelenofTroy.com/emea-en/)