


MEDISANA®


PowerRoll Pro

Art.-No. 79466 | EAN 40 15588 79466 7

- ✓ PowerRoll Soft with depth vibration – also reaches the deeper muscles
- ✓ For activation and regeneration as well as for strength and stability exercises
- ✓ 4 intensity levels
- ✓ LED indication of the selected intensity level



Activation • Regeneration • Fascial Training • Massage



PowerRoll Pro



Art.-No. 79466 | EAN 40 15588 79466 7

Rechargeable Lithium-Ion battery

One-touch operation

4 intensity levels

Surface structure supports the massage action



Application | Features

PowerRoll Pro with depth vibration – also reaches the deeper muscles

For activation and regeneration as well as for strength and stability exercises

4 intensity levels

Simple One-touch operation

LED indication of the selected intensity level

Surface structure supports the massage action

Surface is resistant to moisture

Integrated rechargeable Lithium-Ion battery (charger included)

Max. operating time with fully charged battery: ca. 180 min. (depending on the selected intensity level)

Auto-Stop after 10 minutes of operation

Includes training DVD and training plan

Technical Data

Length: 31 cm, ø15 cm

Weight: approx. 0,75 kg

Battery: Lithium-Ion battery

Package Content

PowerRoll Pro

Training DVD

Training plan

Instruction manual

Charger

Logistical Data

Sales Packaging

1 piece | 15,5 x 36,5 x 15,5 cm | 1,2 kg

Packaging unit

4 pieces | 32,8 x 32,8 x 37,5 cm | 5,5 kg