

therapy



UNITE DE
COMMANDE SUR LE
PRODUIT
WITH THE
OPERATING ELEMENT
ON THE PRODUCT



germanengineering

MG 148 Ceinture de massage Shiatsu

Utilisation versatile et sélective

Massage shiatsu apaisant

Utilisation ciblée pour la nuque, le dos et les épaules

Utilisation facile avec l'unité de commande sur le produit

Choix de 2 directions de massage

Poignée pratique

Utilisation sans les mains possible en attachant la ceinture au corps

Extension de la ceinture incluse

Fonctions lumière et chauffage

Arrêt automatique

Taille : 124 x 14 x 10 cm

Poids : 1,2 kg

Env. 12 watts

Garantie 3 ans

Colisage : 4

Carton de transport : -

EAN-Nr.: 42 11125 64307 2

Art.Nr.: 643.07

MG 148 Shiatsu massage belt

Versatile and selective application

Soothing shiatsu massage

Targeted application for the neck, back and shoulder area

Easy application with operating element on the product

Choice of 2 massage directions

Practical handle

Hands-free application possible by attaching to body

Incl. belt extension

Light and heat function

Automatic switch-off

Measurements (L/H/W): 124 x 14 x 10 cm

Weight: 1.2 kg

Rated output approx. 12 W

3 year guarantee

Sales unit: 4

Shipping carton: -

EAN no.: 42 11125 64307 2

Item number: 643.07

beurer